



TiterZyme ®EIA

human IL-10

Enzyme Immunometric Assay Kit

Catalog No. 900-036

96 Well Kit

Table of Contents

Description	Page	2
Introduction		2
Precautions		2
Materials Supplied		3
Storage		3
Materials Needed but Not Supplied		3
Sample Handling		4
Procedural Notes		5
Reagent Preparation		5
Assay Procedure		6
Calculation of Results		7
Typical Results		7
Typical Standard Curve		8
Performance Characteristics		9
Sample Dilution Recommendations		11
References		11
Limited Warranty		12

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

Description

Assay Designs' human IL-10 TiterZyme® Enzyme Immunometric Assay (EIA) kit is a complete kit for the quantitative determintation of human IL-10 in biological fluids. Please read the complete kit insert before performing this assay. The kit uses a monoclonal antibody to human IL-10 immobilized on a microtiter plate to bind the human IL-10 in the standards or sample. A recombinant human IL-10 Standard is provided in the kit. After a short incubation the excess sample or standard is washed out and a biotinylated monoclonal antibody to human IL-10 is added. This antibody binds to the human IL-10 captured on the plate. After a short incubation the excess antibody is washed out and Streptavidin conjugated to Horseradish peroxidase is added, which binds to the biotinylated human IL-10 antibody. Excess conjugate is washed out and substrate is added. After a short incubation, the enzyme reaction is stopped and the color generated is read at 450 nm. The measured optical density is directly proportional to the concentration of human IL-10 in either standards or samples. For further explanation of the principles and practices of immunoassays please see the excellent books by Chard¹ or Tijssen².

Introduction

Interleukin-10 (IL-10) was initially identified as the 'Cytokine Synthesis Inhibitory Factor' based on the observation that mouse Th1 cytokine production was suppressed by a Th2-derived factor.³ IL-10 is a macrophage deactivating factor, acting on macrophage-monocyte accessory cells to produce its inhibitory effects on T cells and natural killer cells. It also regulates growth and/or differentiation of B cells, mast cells, granulocytes, dendritic cells, karatinocytes and endothelial cells.⁴ The ability of IL-10 to suppress Th1 activities while stimulating Th2 and humoral immune responses has been circumvented by a variety intracellular pathogens known to target macrophages. A number of parasites, bacteria, fungi and viruses depress host immune responses by either inducing host IL-10 production or encode their own IL-10 homologs.⁵ Besides being a potent immunosuppressant, IL-10 is also an antipyretic and functions in mast cell homeostasis.⁶? Circulating levels of IL-10 are increased in allergic asthma, systemic sclerosis, a variety of cancers, post-transplantation patients, and in sepsis.⁸⁻¹² The therapeutic potential of IL-10 includes rheumatoid arthritis, systemic lupus erythematosus, multiple sclerosis, and HIV infections.¹³⁻¹⁷

Precautions

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

- 1. Stop Solution 2 is a 1 normal (1N) hydrochloric acid solution. This solution is caustic; care should be taken in use.
- 2. The activity of the Horseradish peroxidase conjugate is affected by nucleophiles such as azide, cyanide and hydroxylamine.
- 3. We test this kit's performance with a variety of samples, however it is possible that high levels of interfering substances may cause variation in assay results.
- 4. The human IL-10 Standard provided, Catalog No. 80-0658, should be handled with care because of the known and unknown effects of human IL-10.
- 5. The human IL-10 Standard should be stored at or below -20 °C. Do not repeatedly freeze-thaw.

Materials Supplied

- 1. human IL-10 Microtiter Plate, One Plate of 96 Wells, Catalog No. 80-0369

 A plate using break-apart strips coated with monoclonal antibody specific to human IL-10.
- 2. human IL-10 EIA Antibody, 5 mL, Catalog No. 80-0371
 A yellow solution of biotinylated monoclonal antibody to human IL-10.
- 3. Assay Buffer 13, 55 mL, Catalog No. 80-1500 Tris buffered saline containing proteins and detergents.
- 4. human IL-10 EIA Conjugate, 5 mL, Catalog No. 80-0889

 A blue solution of Streptavidin conjugated to Horseradish peroxidase.
- 5. Wash Buffer Concentrate, 100 mL, Catalog No. 80-1287 Tris buffered saline containing detergents.
- 6. human IL-10 Standard, 2 each, Catalog No. 80-0658
 Two vials containing 2,000 pg each of lyophilized recombinant human IL-10.
- 7. TMB Substrate, 10 mL, Catalog No. 80-0350
 A solution of 3,3',5,5' tetramethylbenzidine (TMB) and hydrogen peroxide. Protect from prolonged exposure to light.
- 8. Stop Solution 2, 10 mL, Catalog No. 80-0377
 A 1N solution of hydrochloric acid in water. Keep tightly capped. Caution: Caustic.
- 9. human IL-10 Assay Layout Sheet, 1 each, Catalog No. 30-0105
- 10. Plate Sealer, 3 each, Catalog No. 30-0012

Storage

All components of this kit, **except the Standards**, are stable at 4 °C until the kit's expiration date. The Standards **must** be stored at or below -20 °C.

Materials Needed but Not Supplied

- 1. Deionized or distilled water.
- 2. Precision pipets for volumes between 100 μ L and 1,000 μ L.
- 3. Repeater pipet for dispensing 100 μL.
- 4. Disposable beakers for diluting buffer concentrates.
- 5. Graduated cylinders.
- 6. A microplate shaker.
- 7. Adsorbent paper for blotting.
- 8. Microplate reader capable of reading at 450 nm, preferably with correction between 570 nm and 590 nm.
- 9. Graph paper for plotting the standard curve.

Sample Handling

Assay Designs' TiterZyme® EIA is compatible with human IL-10 samples in a wide range of matrices. Samples diluted sufficiently into the proper diluent can be read directly from a standard curve. Please refer to the Sample Recovery recommendations on page 11 for details of suggested dilutions.

Culture fluids, serum and plasma are suitable for use in the assay. Samples containing a visible precipitate must be clarified prior to use in the assay. Do not use grossly hemolyzed or lipemic specimens. Samples in the majority of tissue culture media, including those containing fetal bovine serum, can also be read in the assay, provided the standards have been diluted into the tissue culture media instead of Assay Buffer 13. There will be a small change in binding associated with running the standards and samples in media. Users should only use standard curves generated in media or buffer to calculate concentrations of human IL-10 in the appropriate matrix.

Samples must be stored frozen to avoid loss of bioactive human IL-10. If samples are to be run within 24 hours, they may be stored at 4 °C. Otherwise, samples must be stored frozen at or below -70 °C to avoid loss of bioactive human IL-10. Excessive freeze/thaw cycles should be avoided. Prior to assay, frozen sera should be brought to room temperature slowly and gently mixed by hand. Do not thaw samples in a 37 °C incubator. Do not vortex or sharply agitate samples.

Procedural Notes

- 1. Do not mix components from different kit lots or use reagents beyond the kit expiration date.
- 2. Allow all reagents to warm to room temperature for at least 30 minutes before opening.
- 3. Standards can be made up in either glass or plastic tubes.
- 4. Pre-rinse the pipet tip with reagent, use fresh pipet tips for each sample, standard and reagent.
- 5. Pipet standards and samples to the bottom of the wells.
- 6. Add the reagents to the side of the well to avoid contamination.
- 7. This kit uses break-apart microtiter strips, which allow the user to measure as many samples as desired. Unused wells must be kept desiccated at 4 °C in the sealed bag provided. The wells should be used in the frame provided.
- 8. Prior to addition of antibody, conjugate and substrate, ensure that there is no residual wash buffer in the wells. Any remaining wash buffer may cause variation in assay results.
- 9. It is important that the matrix for the standards and samples be as similar as possible. Human IL-10 samples diluted with Assay Buffer 13 should be run with a standard curve diluted in the same buffer. Serum and plasma samples and samples with a protein concentration greater than 10% should be evaluated against a standard curve run in Assay Buffer 13 while tissue culture samples should be read against a standard curve diluted in the same complete but non-conditioned media. See Reagent Preperation, step #2.

Reagent Preparation

1. Wash Buffer

Prepare the Wash Buffer by diluting 50 mL of the supplied concentrate with 950 mL of deionized water. This can be stored at room temperature until the kit expiration, or for 3 months, whichever is earlier.

2. human IL-10 Standards

Reconstitute one vial of human IL-10 Standard with 2 mL standard diluent (Assay Buffer 13 or Tissue Culture Media). Mix thoroughly without foaming. Label seven 12x75 mm glass tubes #1 through #7. Pipet 500 μ L of standard diluent into tubes #1 through #7. Add 500 μ L of the 1,000 pg/mL Standard to tube #1. Vortex thoroughly. Add 500 μ L of tube #2 and vortex thoroughly. Add 500 μ L of tube #2 to #3 and vortex thoroughly. Continue this for tubes #4 through #7.

The concentration of human IL-10 in tubes #1 through #7 will be 500, 250, 125, 62.5, 31.25, 15.62 and 7.81 pg/mL, respectively. See human IL-10 Assay Layout Sheet for dilution details. STORE STANDARD AT -20 °C, avoid repeated freeze-thaws.

Assay Procedure

Bring all reagents to room temperature for at least 30 minutes prior to opening.

All standards, controls and samples should be run in duplicate.

- 1. Refer to the Assay Layout Sheet to determine the number of wells to be used and put any remaining wells with the desiccant back into the pouch and seal the ziploc. Store unused wells at 4 °C.
- 2. Pipet $100 \,\mu\text{L}$ of standard diluent (Assay Buffer 13 or Tissue Culture Media) into the S0 (0 pg/mL standard) wells.
- 3. Pipet 100 µL of Standards #1 through #7 into the appropriate wells.
- 4. Pipet 100 μL of the Samples into the appropriate wells.
- 5. Tap the plate gently to mix the contents.
- 6. Seal the plate and incubate at room temperature on a plate shaker for 1 hour at ~500 rpm.
- 7. Empty the contents of the wells and wash by adding 400 µL of wash solution to every well. Repeat the wash 3 more times for a total of **4 washes**. After the final wash, empty or aspirate the wells and firmly tap the plate on a lint free paper towel to remove any remaining wash buffer.
- 8. Pipet 100 µL of yellow Antibody into each well, except the Blank.
- 9. Seal the plate and incubate at room temperature on a plate shaker for 1 hour at ~500 rpm.
- 10. Empty the contents of the wells and wash by adding 400 μL of wash solution to every well. Repeat the wash 3 more times for a total of **4 washes**. After the final wash, empty or aspirate the wells and firmly tap the plate on a lint free paper towel to remove any remaining wash buffer.
- 11. Add 100 µL of blue Conjugate to each well, except the Blank.
- 12. Seal the plate and incubate at room temperature on a plate shaker for 30 minutes at ~500 rpm.
- 13. Empty the contents of the wells and wash by adding 400 µL of wash solution to every well. Repeat the wash 3 more times for a total of **4 washes**. After the final wash, empty or aspirate the wells and firmly tap the plate on a lint free paper towel to remove any remaining wash buffer.
- 14. Pipet 100 μL of Substrate Solution into each well.
- 15. Incubate for 15 minutes at room temperature on a plate shaker at ~500 rpm.
- 16. Pipet $100 \mu L$ Stop Solution 2 to each well. This stops the reaction and the plate should be read immediately.
- 17. Blank the plate reader against the Blank wells, read the optical density at 450 nm, preferably with correction between 570 and 590 nm. If the plate reader is not able to be blanked against the Blank wells, manually subtract the mean optical density of the Blank wells from all the readings.

Calculation of Results

Several options are available for the calculation of the concentration of human IL-10 in the samples. We recommend that the data be handled by an immunoassay software package utilizing a 4 parameter logistic curve fitting program. If data reduction software is not readily available, the concentration of human IL-10 can be calculated as follows:

- 1. Calculate the average net Optical Density (OD) bound for each standard and sample by subtracting the average Blank OD from the average OD for each standard and sample.
 - Average Net OD = Average OD Average Blank OD
- 2. Using linear graph paper, plot the Average Net OD for each standard versus human IL-10 concentration in each standard. Approximate a straight line through the points. The concentration of human IL-10 in the unknowns can be determined by interpolation.


Typical Results

The results shown below are for illustration only and **should not** be used to calculate results from another assay.

Sample	Average OD	Net OD	h IL-10 <u>(pg/mL)</u>
Blank	0.073		
S0	0.087	0.014	0
S 1	1.978	1.905	500
S2	0.981	0.908	250
S3	0.528	0.455	125
S4	0.302	0.229	62.5
S5	0.193	0.120	31.25
S6	0.139	0.066	15.62
S7	0.103	0.030	7.81
Unknown 1	1.699	1.626	430.0
Unknown 2	0.427	0.354	97.2

Typical Standard Curve

Typical standard curves are shown below. These curves **must not** be used to calculate human IL-10 concentrations; each user must run a standard curve for each assay.



human IL-10 Conc. (pg/mL)

Performance Characteristics

The following parameters for this kit were determined using the guidelines listed in the National Committee for Clinical Laboratory Standards (NCCLS) Evaluation Protocols¹⁸.

Sensitivity

Sensitivity was calculated by determining the average optical density bound for sixteen (16) wells run with 0 pg/mL Standard, and comparing to the average optical density for sixteen (16) wells run with Standard #7. The detection limit was determined as the concentration of human IL-10 measured at two (2) standard deviations from the 0 pg/mL Standard along the standard curve.

Mean OD for S0 = Mean OD for Standard #7 =	$0.009 \pm 0.006 (14.61\%)$ $0.034 \pm 0.005 (13.21\%)$
Delta Optical Density $(7.81 - 0 \text{ pg/mL}) = 0.034$	4 - 0.009 = 0.025
2 SD's of 0 pg/mL Standard = 2 x 0.006 =	0.012
Sensitivity = $\frac{0.012}{0.025}$ x $7.81 \text{ pg/mL} = \frac{0.012}{0.025}$	3.75 pg/mL

Linearity

A sample containing 297.4 pg/mL human IL-10 was serially diluted 5 times 1:2 in the Assay Buffer 13 supplied in the kit and measured in the assay. The data was plotted graphically as actual human IL-10 concentration versus measured human IL-10 concentration.

The line obtained had a slope of 1.01 with a correlation coefficient of 0.998.

Precision

Intra-assay precision was determined by taking samples containing low, medium and high concentrations of human IL-10 and running these samples multiple times (n=16) in the same assay. Interassay precision was determined by measuring three samples with low, medium and high concentrations of human IL-10 in multiple assays (n=8).

The precision numbers listed below represent the percent coefficient of variation for the concentrations of human IL-10 determined in these assays as calculated by a 4 parameter logistic curve fitting program.

	h IL-10	Intra-assay	Inter-assay
	(pg/mL)	<u>% CV</u>	<u>% CV</u>
	100 5	~ .	
Low	102.7	5.4	
Medium	194.3	3.9	
High	435.9	2.2	
Low	106.1		4.7
Medium	203.2		4.5
High	459.4		3.7

Cross Reactivities

The TiterZyme® human IL-10 EIA kit is specific for bioactive human IL-10. It is unaffected by the presence of recombinant mouse IL-10 and the following recombinant human proteins: IL-1α, IL-1β, IL-2, IL-3, IL-4, IL-6, IL-7, IL-8, IL-11, TNF-α, TNF-β, GM-CSF and SCF.

Sample Recoveries

Please refer to pages 4 and 5 for Sample Handling recommendations and Standard preparation.

Human IL-10 concentrations were measured in human plasma and serum, and tissue culture media. Human IL-10 was spiked into the undiluted samples of these matrices which were then diluted with the appropriate diluent and assayed in the kit. The following results were obtained:

		Recommended
Sample	% Recovery*	Dilution *
Human heparin Plasma	103.0	None
Human Serum	86.7	None
Tissue Culture Media	95.2	None

^{*} See Sample Handling instructions on page 4 for details.

References

- 1. T. Chard, "An Introduction to Radioimmunoassay & Related Techniques, 4th Ed.", (1990) Amsterdam: Elsevier.
- 2. P. Tijssen, "Practice & Theory of Enzyme Immunoassays", (1985) Amsterdam: Elsevier.
- 3. D.F. Fiorentino, et al., J. Exp. Med., (1989) 170: 2081.
- 4. K.W. Morre, et al., Annu. Rev. Immunol. (2001) 19: 683.
- 5. S. Redpath, et al., <u>Trends Microbiol.</u>, (2001) <u>9/2</u>: 86.
- 6. W. Kozak, et al., Ann. NY Acad. Sci., (2000), 917: 121.
- 7. C.P. Shelburne & J.J. Ryan., <u>Immunol. Rev.</u>, (2001) <u>179</u>: 82.
- 8. C.K. Wong, et al., Clin. Exp. I mmunol., (2001) 125/2: 177.
- 9. S. Sato, et al., J. Dermatol. Sci., (2001) 27/2: 140.
- 10. S. Mocellin, et al., J Immunother., (2001) 24/5: 392.
- 11. L. Kubala, et al., Cytokine., (2001) 16/3: 97.
- 12. M.R. Pinsky, Contrib. Nephrol., (2001) 132: 354.
- 13. A. Kumar & W.D. Cheery, <u>Arch. Immunol. Ther. Exp.</u>, (2000) <u>48/6</u>: 529.
- 14. B. Joe, et al., Curr. Rheumatol. Rep., (1999) 1/2: 139.
- 15. F.C. Arnett, Int. Rev. Immunol., (2000) 19/4-5: 297.
- 16. C.L. Carp, et al., Curr. Opin. Neurol., (2000) <u>14/3</u>; 361.
- 17. M.A. Ostrowski, et al., J. Infect. Dis., (2001) 184/10: 1268.
- 18. National Committee for Clinical Laboratory Standards Evaluation Protocols, SC1, (1989) Villanova, PA: NCCLS.

LIMITED WARRANTY

Assay Designs, Inc. warrants that at the time of shipment this product is free from defects in materials and workmanship. This warranty is in lieu of any other warranty expressed or implied, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose.

Assay Designs must be notified of any breach of this warranty within 48 hours of receipt of the product. No claim shall be honored if Assay Designs is not notified within this time period, or if the product has been stored in any way other than outlined in this publication. The sole and exclusive remedy of the customer for any liability based upon this warranty is limited to the replacement of the product, or refund of the invoice price of the goods.



For more details concerning the information within this kit insert, or to order any of Assay Designs' products, please call (734) 668-6113 between 8:30 a.m. and 5:30 p.m. EST. Orders or technical questions can also be transmitted by fax or e-mail 24 hours a day.

Material Safety Data Sheet (MSDS) available on our website or by fax.

Assay Designs, Inc. Telephone: (734) 668-6113

5777 Hines Drive (800) 833-8651 (USA & Canada only)

Ann Arbor, MI 48108 Fax: (734) 668-2793

U.S.A. e-mail: info@assaydesigns.com website: www.assaydesigns.com

Simplify Your Science®

Catalog No. 25-0200

© 2000



April 26, 2007

