



TiterZyme[®] EIA

human Cyclooxygenase-II Enzyme Immunometric Assay Kit

Catalog No. 900-094

96 Determination Kit

Table of Contents

Description	Page	2
Introduction		2
Precautions		2
Materials Supplied		3
Storage		3
Materials Needed but Not Supplied		3
Sample Handling		4
Procedural Notes		4
Reagent Preparation		5
Assay Procedure		6
Calculation of Results		7
Typical Results		7
Typical Standard Curve		8
Performance Characteristics		9
Sample Dilution Recommendations		11
References		11
Limited Warranty		12


FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

Description

Assay Designs' human COX-II Enzyme Immunometric Assay (EIA) kit is a complete kit for the quantitative determination of human COX-II in cell culture lysates. Please read the complete kit insert before performing this assay. The kit uses a monoclonal antibody to human COX-II immobilized on a microtiter plate to bind the human COX-II in the standard or sample. A recombinant human COX-II Standard is provided in the kit. After a short incubation the excess standard or sample is washed out and a rabbit polyclonal antibody to human COX-II labeled with the enzyme Horseradish peroxidase is added. This labeled antibody binds to the human COX-II captured on the plate. After a short incubation the excess labeled antibody is washed out and substrate is added. The substrate reacts with the labeled antibody bound to the human COX-II captured on the plate. After a short incubation, the enzyme reaction is stopped and the color generated is read at 450 nm. The measured optical density is directly proportional to the concentration of human COX-II in either standards or samples. For further explanation of the principles and practice of immunoassays please see the excellent books by Chard¹ or Tijssen².

Introduction

Cyclooxygenase (COX, also known as Prostaglandin G/H synthase) is a membrane bound enzyme responsible for the oxidation of arachidonic acid to Prostaglandin G₂ (PGG₂) and the subsequent reduction of PGG₂ to PGH₂^{3,4}. The conversion is shown below. These reactions are the first steps in the formation of a variety of prostanoids. COX has been shown to be expressed in at least two different isoforms, a constitutively expressed form, COX-I, and an inducible form, COX-II. COX-I is thought to regulate a number of 'housekeeping' functions, such as vascular hemostasis, renal blood flow, and maintenance of glomerular function⁵. Inflammation mediators such as growth factors, cytokines and endotoxin induce COX-II expression in a number of cellular systems^{6,7}. The effect of various non-steroidal anti-inflammatory drugs (NSAID's) on the activity of COX-I and -II is an area of considerable interest⁸.



Precautions

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

1. Stop Solution is a 1 normal (1N) sulfuric acid solution. This solution is caustic; care should be taken in use.
2. The activity of the Horseradish peroxidase conjugate is affected by nucleophiles, such as azide, cyanide and hydroxylamine.
3. We test this kit's performance with a variety of samples, however, it is possible that high levels of interfering substances may cause variation in assay results.
4. The human COX-II Standard provided, Catalog No. 80-0908, should be handled with care because of the known and unknown effects of cyclooxygenase.
5. The human COX-II Standard and Labeled Antibody should be stored at -20°C. Do not repeatedly freeze-thaw.

Materials Supplied

1. **human COX-II Microtiter Plate, One Plate of 96 Wells, Catalog No. 80-0906**
A strip microtiter plate coated with mouse antibody specific to human COX-II.
2. **human COX-II Labeled Antibody, 1 vial, Catalog No. 80-0907**
Rabbit antibody to human COX-II conjugated to Horseradish peroxidase.
3. **Assay Buffer, 30 mL, Catalog No. 80-0170**
Phosphate buffered saline containing proteins and detergents.
4. **Labeled Antibody Diluent, 10.5 mL, Catalog No. 80-0182**
Phosphate buffered saline containing proteins and detergents.
5. **Wash Buffer Concentrate, 50 mL, Catalog No. 80-0171**
Phosphate buffered saline containing detergents.
6. **human COX-II Standard, 1 vial, Catalog No. 80-0908**
A vial containing 275 ng of recombinant human COX-II.
7. **TMB Substrate, 15 mL, Catalog No. 80-1342**
A solution of 3,3',5,5' tetramethyl benzidine (TMB) and hydrogen peroxide. Ready to use.
8. **Stop Solution, 11 mL, Catalog No. 80-0176**
A 1N solution of sulfuric acid in water. Keep tightly capped. Caution: **Caustic**.
9. **human COX-II Assay Layout Sheet, 1 each, Catalog No. 30-0160**
10. **Plate Sealer, 2 each, Catalog No. 30-0012**

Storage

All components of this kit, are stable at 4°C until the kit's expiration date.

Materials Needed but Not Supplied

1. Deionized or distilled water. No difference in assay results are seen with distilled water.
2. Precision pipets for volumes between 100 μ L and 1,000 μ L.
3. Repeater pipet for dispensing 100 μ L.
4. Disposable beaker for diluting buffer concentrates.
5. Graduated cylinders.
6. A 37°C incubator.
7. Adsorbent paper for blotting.
8. Microplate reader capable of reading at 450 nm, preferably with correction between 570 and 590 nm.
9. Graph paper for plotting the standard curve.

Sample Handling

Assay Designs' TiterZyme® EIA is compatible with human COX-II samples in cell culture lysates only. This assay is not for use in serum or plasma. Samples diluted sufficiently into Assay Buffer ($\geq 1:2$) can be read directly from the standard curve. Please refer to the Sample Recovery recommendations on page 11 for details of suggested dilutions.

Samples can be prepared by sonicating cells in TNE (10 mM Tris, pH 8.0, 0.15 M NaCl, 1% NP-40, 1mM EDTA) or RIPA (25 mM Tris, pH 7.4, 0.15 M KCl, 1% NP-40, 5 mM EDTA, 0.5% Sodium deoxycholate, 0.1% SDS) buffer. Five cycles of 30 second bursts at 1 minute intervals on ice is usually adequate for disruption. For optimal results, a small aliquot of the lysate should be examined under a microscope to verify sample disruption. The supernatant is collected after a 5 minute centrifugation at 15,000 rpm, 4°C and assayed immediately. Typically, the supernatant from a lysate of 10^7 cells in 0.5 mL buffer will need an additional 5 to 10 fold dilution prior to the assay. Buffers containing CHAPS and sucrose monolaurate should be avoided.

Plasma and serum are not suitable matrices for the measurement of human COX-II due to their low circulating concentrations.

Procedural Notes

1. Do not mix components from different kit lots or use reagents beyond the kit expiration date.
2. Allow all reagents to warm to room temperature for at least 30 minutes before opening.
3. Standards can be made up in either glass or plastic tubes.
4. Pre-rinse the pipet tip with reagent, use fresh pipet tips for each sample, standard and reagent.
5. Pipet standards and samples to the bottom of the wells.
6. Add the reagents to the side of the well to avoid contamination.
7. This kit uses plates with removable strips. Unused strips must be kept desiccated at 4°C in the sealed foil bag. The strips should be used in the frame provided.
8. **Prior to addition of standard, antibody, and substrate, ensure that there is no residual wash buffer in these wells. Any remaining wash buffer may cause variation in assay results.**

Reagent Preparation

1. Wash Buffer

Prepare Wash Buffer by diluting 25 mLs of the supplied concentrate with 975 mLs of deionized water. This can be stored at 4°C until the kit expiration, or for 3 months, whichever is earlier.

2. human COX-II Standard

Add 500 µL of deionized water to the human COX-II Standard. Let it sit at room temperature for 5 minutes. Mix it gently. This solution contains 550 ng/mL human COX-II.

Label seven 12 x 75 mm glass tubes #1 through #7. Pipet 750 µL of Assay Buffer into tube #1. Pipet 500 µL of Assay Buffer into tubes #2 through #7. Add 250 µL of the 550 ng/mL standard to tube #1. Vortex. Add 500 µL of tube #1 to tube #2 and vortex thoroughly. Continue this for tubes #3 through #7.

The concentration of human COX-II in tubes #1 through #7 will be 137.5, 68.75, 34.38, 17.19, 8.59, 4.3 and 2.15 ng/mL respectively. See human COX-II Assay Layout Sheet for dilution details. STORE STANDARD AT -20°C, avoid repeated freeze/thaws.

3. Preparation of Labeled Antibody Conjugate

Add the entire contents of one (1) bottle of Labeled Antibody Diluent to the vial of human COX-II Antibody Conjugate. Let it stand at room temperature for 5 minutes and then vortex it gently. After reconstitution any unused Labeled Antibody should be aliquoted and stored at -20°C.

Assay Procedure

Bring all reagents to room temperature for at least 30 minutes prior to opening.

All standards and samples should be run in duplicate.

1. Refer to the Assay Layout Sheet to determine the number of wells to be used and put any remaining strips with the desiccant back into the foil pouch and seal the ziploc. Store unused strips at 4°C.
2. Pipet 100 µL of Assay Buffer into the S0 (0 ng/mL Standard) wells.
3. Pipet 100 µL of Standards #1 through #7 into the appropriate wells.
4. Pipet 100 µL of the Samples into the appropriate wells.
5. Tap the plate gently to mix the contents.
6. Seal the plate and incubate at 37°C for 1 hour.
7. Empty the contents of the wells and wash by adding 200 µL of wash solution to every well. Repeat the wash 6 more times for a total of **7 washes**. After the final wash, empty or aspirate the wells, and firmly tap the plate on a lint-free paper towel to remove any remaining wash buffer.
8. Pipet 100 µL of the Labeled Antibody into each well, except the Blank.
9. Seal the plate and incubate at 4°C for 30 minutes. Prepare Substrate (See page 5, Section 5).
10. Empty the contents of the wells and wash by adding 200 µL of wash solution to every well. Repeat the wash 8 more times for a total of **9 washes**. After the final wash, empty or aspirate the wells, and firmly tap the plate on a lint-free paper towel to remove any remaining wash buffer.
11. Add 100 µL of the TMB Substrate to each well.
12. Incubate for 30 minutes at room temperature in the dark.
13. Add 100 µL of Stop Solution to each well.
14. Blank the plate reader against the Blank wells, read the optical density at 450nm, preferably with correction between 570 and 590 nm. If the plate reader is not able to be blanked against the Blank wells, manually subtract the mean optical density of the blank wells from all readings.

Calculation of Results

Several options are available for the calculation of the concentration of human COX-II in the samples. We recommend that the data be handled by an immunoassay software package utilizing a 4 parameter logistic curve fitting program. If data reduction software is not readily available, the concentration of human COX-II can be calculated as follows:

1. Calculate the average net Optical Density (OD) bound for each standard and sample by subtracting the average Blank OD from the average OD for each standard and sample.

$$\text{Average Net OD} = \text{Average OD} - \text{Average Blank OD}$$

2. Using linear graph paper, plot the Average Net OD for each standard versus human COX-II concentration in each standard. Approximate a straight line through the points. The concentration of human COX-II in the unknowns can be determined by interpolation.


Typical Results

The results shown below are for illustration only and **should not** be used to calculate results from another assay.

<u>Sample</u>	<u>Average OD</u>	<u>Net OD</u>	human COX-II (ng/mL)
Blank	0.050		
S0	0.116	0.066	0
S1	2.420	2.370	137.5
S2	1.765	1.715	68.75
S3	1.137	1.087	34.38
S4	0.659	0.609	17.19
S5	0.419	0.369	8.59
S6	0.260	0.210	4.30
S7	0.185	0.135	2.15
Unknown 1	2.069	2.019	104.45
Unknown 2	0.877	0.827	32.92

Typical Standard Curve

The typical standard curve shown below **must not** be used to calculate human COX-II concentrations; the user must run a standard curve for each plate used.



Performance Characteristics

The following parameters for this kit were determined using the guidelines listed in the National Committee for Clinical Laboratory Standards (NCCLS) Evaluation Protocols⁹.

Sensitivity

Sensitivity was calculated by determining the average optical density bound for sixteen (16) wells run with 0 ng/mL Standard and comparing to the average optical density for sixteen (16) wells run with Standard #7. The detection limit was determined as the concentration of human COX-II measured at two (2) standard deviations from the 0 ng/mL Standard along the standard curve.

Average Optical Density for the S0 = 0.067 ± 0.004 (6.5%)
Average Optical Density for Standard #7 = 0.136 ± 0.003 (1.9%)

Delta Optical Density (2.15-0 ng/mL)= 0.136 - 0.067 = 0.069

2 SD's of the 0 ng/mL Standard = 2 x 0.004 = 0.008

Sensitivity = $\frac{0.008}{0.069} \times 2.15 \text{ ng/mL} = \mathbf{0.249 \text{ ng/mL}}$

Linearity

A sample containing 104.45 ng/mL human COX-II was diluted 5 times 1:2 into the kit Assay Buffer and measured in the assay. The data was plotted graphically as actual human COX-II concentration versus measured human COX-II concentration.

The line obtained had a slope of 0.8955 and a correlation coefficient of 1.000.

Precision

Intra-assay precision was determined by taking samples containing low, medium and high concentrations of human COX-II and running these samples multiple times (n=4) in the same assay. Inter-assay precision was determined by measuring three samples with low, medium and high concentrations of human COX-II in multiple assays (n=7).

The precision numbers listed below represent the percent coefficient of variation for the concentrations of human COX-II determined in these assays as calculated by a 4 parameter logistic curve fitting program.

	<u>human COX-II</u> <u>(ng/mL)</u>	<u>Intra-assay</u> <u>%CV</u>	<u>Inter-assay</u> <u>%CV</u>
Low	6.4	7.6	
Medium	52.8	5.6	
High	139.9	2.5	
Low	6.2		4.8
Medium	52.5		5.9
High	137.9		6.0

Cross Reactivities

The human COX-II TiterZyme® EIA kit is specific for human COX-II. The cross reactivity to human COX-I is <0.1%.

Sample Recoveries

Please refer to pages 4 and 5 for Sample Handling recommendations and Standard preparation.

Human COX-II concentrations were measured in TNE buffer. Human COX-II was spiked into the undiluted sample which was then diluted with the kit Assay Buffer and assayed in the kit. The following results were obtained:

<u>Sample</u>	<u>% Recovery</u> *	<u>Recommended Dilution</u> *
TNE Buffer	103.8	≥1:2

* See Sample Handling instructions on page 4 for details.

References

1. T. Chard, in "An Intro. to Radioimmunoassay & Related Tech.", (1990), 4th Ed., Elsevier, Amsterdam.
2. P. Tijssen, in "Practice and Theory of Enz. Immunoassays", (1985), Elsevier, Amsterdam.
3. J.R. Vane and R.M. Botting, Inflamm. Res., (1995), 44, 1-10.
4. W.L. Smith, L.J. Marnett and D.L. DeWitt in "Intracellular Messengers", (1993), 101-136, Pergamon Press, London, C.W. Taylor, (Ed).
5. P. Mene, M.S. Simonson, and M.J. Dunn, Physiol. Rev., (1989), 69, 1347-1423.
6. H. Herschman, B. Fletcher and D. Kujubu, J. Lipid. Med., (1993), 6, 89-99.
7. E.A. Meade, W.L. Smith, and D.L. DeWitt, ibid., (1993), 6, 119-129.
8. S. Kargman, et al., Biochem. Pharm., (1996), 52, 1113-1125.
9. National Committee for Clinical Laboratory Standards Evaluation Protocols, SC1, 1989, NCCLS, Villanova, PA, 19085.

LIMITED WARRANTY

Assay Designs, Inc. warrants that at the time of shipment this product is free from defects in materials and workmanship. This warranty is in lieu of any other warranty expressed or implied, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose.

Assay Designs must be notified of any breach of this warranty within 48 hours of receipt of the product. No claim shall be honored if Assay Designs is not notified within this time period, or if the product has been stored in any way other than outlined in this publication. The sole and exclusive remedy of the customer for any liability based upon this warranty is limited to the replacement of the product, or refund of the invoice price of the goods.



For more details concerning the information within this kit insert, or to order any of Assay Designs' products, please call (734) 668-6113 between 8:30 a.m. and 5:30 p.m. EST. Orders or technical questions can also be transmitted by fax or e-mail 24 hours a day.

Material Safety Data Sheet (MSDS) available on our website or by fax.

**Assay Designs, Inc.
5777 Hines Drive
Ann Arbor, MI 48108
U.S.A.**

**Telephone: (734) 668-6113
(800) 833-8651 (USA & Canada only)
Fax: (734) 668-2793
website: www.assaydesigns.com**

Simplify Your Science®

Catalog No. 25-0409

May 22, 2007

© 2001

