

KT117 - Anti Human Macrophage (CD68) Monoclonal Antibody (Clone No. PM-1K)

Macrophages are present in nearly all tissues and organs of the body. They are differentiated from monocytes derived from the bone marrow. Macrophages and monocytes are phagocytes, acting in both innate immunity and cell-mediated immunity of vertebrate animals. Their function is to phagocytize cellular debris and pathogens, and to stimulate lymphocytes and other immune cells to respond to the pathogen.

Macrophages are able to be identified immunohistochemically by virtue of the presence of monocyte/macrophage-associated antigens such as CD68. Antibodies recognizing CD68 have been used as some of the best reagents to detect macrophages in tissues.



CD68 is a 110-kD transmembrane glycoprotein that is highly expressed by human monocytes and tissue macrophages. It is a member of the lysosomal/endosomal-associated membrane glycoprotein (LAMP) family.

This antibody recognizes CD68 antigen. The molecular size of the antigen identified by this antibody was 110 kDa. Immunoprecipitated antigen by this antibody was also recognized by the other CD68 antibodies such as KP-1 and PG-M1.

In immunohistochemical assays, this antibody recognizes freshly isolated human blood monocytes and tissue macrophages. This antibody also recognizes macrophages obtained from guinea pigs, pigs, bovine species, and monkeys. Since this antibody strongly labels guinea pig macrophages, this antibody will be suitable to examine such macrophages in experimental guinea pig models.

This antibody will be very useful to research of CD68, macrophage, allergic diseases and delayed hypersensitivity.

Package Size	50µg (200µL/ vial)
Format	Mouse monoclonal antibody 0.25mg/mL
Buffer	PBS [containing 2% Block Ace as a stabilizer, 0.1%Proclin as a bacteriostat]
Storage	Store below -20°C Once thawed, store at 4°C. Repeated freeze-thaw cycles should be avoided.
Subclass	IgG2b, κ
Purification method	The splenic lymphocytes from BALB/c mouse, immunized with human peritoneal cells from patients with endometriosis incubated for 24 hours, were fused with mouse NS-1 myeloma cells. The hybridoma cell line with positive reaction was grown in ascitic fluid of BALB/c mouse, from which the antibody was purified by Protein G affinity chromatography.
Working dilution for	immunohistochemistry: 10µg/mL



Left; Human spleen (frozen section): Red pulp macrophages are positive.

Right; Small Intestine of Guinea Pig (paraffin section): Macrophages in lamina propria are positive.

Takeya M., Department of Cell Pathology, Graduate School of Medical Sciences, Kumamoto University, Kumamoto, Japan

Gentaur Europe

Avenue de l'Armée 68
B-1040 Brussels
Tel: +32 16 58 90 45
Fax: +32 16 50 90 45
info@gentaur.com
www.gentaur.com

Gentaur Germany

IME, Forckenbeckstrasse 6
D-52074 Aachen Germany
Tel: 0241 6085 13140
Fax: 0241 6085 33033
info@bioprice.com
www.bioprice.com

Gentaur France

9, Rue Lagrange
75005 Paris, France
Tel: +33(0)143 25 0150
Fax: +33(0)143 25 0160
info@clonagen.com
www.clonagen.com

Gentaur Italy

Tel: 02 36 00 65 93
Fax: 02 36 00 65 94
info@labprice.com
www.labprice.com



GENTAUR

□Distribution of positive reactivities of PM-1K with human monocyte/macrophages□

Tissues	Cells	Tissues	Cells
Heart	Intramuscular Mφ	Lymph nodes	Mφ in follicles Mφ in paracortical areas
Lung	Alveolar Mφ	Pancreas, Salivary, Thyroid, Adrenals, Urinary bladder, Prostate, Trachea	Interstitial Mφ Stomach, Small and Large intestines Mφ in lamia propria
Liver	Kupffer cells Mφ in portal triads	Skin	Dermal Mφ Langerhans cells
Kidney	Interstitial Mφ Uriniferous tubule cells	Blood monocytes	Freshly isolated monocytes
Spleen	Red pulp Mφ White pulp Mφ	Blood neutrophils	Freshly isolated neutrophils
Thymus	Mφ in cortex Mφ in medulla	Myeloid cell lines	THP-1 MonoMac6

Mφ : macrophage

□Reference□*Application Reference

1. Horikawa T, Komohara Y, Kiyota E, Terasaki Y, Takagi K, Takeya M. .Detection of guinea pig macrophages by a new CD68 monoclonal antibody, PM-1K. *J Mol Histol*, 37:15-25, 2006*
2. Suenaga Y, Katabuchi H, Fukumatsu Y, Okamura H. Distribution and cytological properties of macrophages in human Fallopian tubes. *Acta Anat (Basel)*, 163:10-19, 1998
3. Imamura T, Iyama K, Takeya M, Kambara T, Nakamura S. Role of macrophage tissue factor in the development of the delayed hypersensitivity reaction in monkey skin. *Cell Immunol*, 152:614-622, 1993
4. Okamoto M, Yamamoto T, Matsubara S, Kukita I, Takeya M, Miyauchi Y, Kambara T. Factor XIII-dependent generation of 5th complement component(C5)-derived monocyte chemotactic factor coinciding with plasma clotting. *Biochem Biophys Acta*, 1138:53-61, 1992

Gentaur Europe
Avenue de l'Armée 68
B-1040 Brussels
Tel: +32 16 58 90 45
Fax: +32 16 50 90 45
info@gentaur.com
www.gentaur.com

Gentaur Germany
IME, Forckenbeckstrasse 6
D-52074 Aachen Germany
Tel: 0241 6085 13140
Fax: 0241 6085 33033
info@bioprice.com
www.bioprice.com

Gentaur France
9, Rue Lagrange
75005 Paris, France
Tel: +33(0)143 25 0150
Fax: +33(0)143 25 0160
info@clonagen.com
www.clonagen.com

Gentaur Italy
Tel: 02 36 00 65 93
Fax: 02 36 00 65 94
info@labprice.com
www.labprice.com