



[www.clonagen.com](http://www.clonagen.com)

[www.bioxys.com](http://www.bioxys.com)

[www.gentaur.com](http://www.gentaur.com)

## QuantiChrom™ ATPase/GTPase Assay Kit (DATG-200)

### DESCRIPTION

ATPases and GTPases catalyze the decomposition of ATP or GTP into ADP or GDP and free phosphate ion. These enzymes play key roles in transport, signal transduction, protein biosynthesis and cell differentiation. BioAssay Systems' QuantiChrom™ ATPase/GTPase Assay Kit offers a highly sensitive method for determining ATPase/GTPase activities in a microplate format. Its proprietary formulation features a single reagent for accurate determination of enzyme activity in 30 min at room temperature. The improved malachite green reagent forms a stable dark green color with liberated phosphate, which is measured on a plate reader (600 - 660 nm).

### KEY FEATURES

**High sensitivity:** detection of as little of 2 pmoles of free phosphate.

**Fast and convenient:** single reagent, homogeneous "mix-and-measure" assay allows quantitation of enzyme activity within 30 minutes.

**Robust and amenable to HTS:** detection at 620nm greatly reduces potential interference by colored compounds. Z' factors of >0.7 are observed in 96-well and 384-well plates. Can be readily automated on HTS liquid handling systems.

### APPLICATIONS

**Determination** of ATPase and GTPase activity.

**Drug Discovery:** high-throughput screen for ATPase or GTPase inhibitors.

### KIT CONTENTS: 200 ASSAYS IN 96-WELL PLATE

**Reagent:** 50 mL **Assay Buffer:** 10 mL

**Standard:** 1mL 1 mM phosphate

**Storage conditions.** The reagents and standard are stable for one year when stored at 4°C.

**Precautions:** reagent contains 0.27 M H<sub>2</sub>SO<sub>4</sub>. Normal precautions for laboratory reagents should be exercised while using the reagents. Please refer to Material Safety Data Sheet for detailed information.

**Important:** All reagents must be brought to room temperature before use. Before each assay, it is important to check that enzyme preparations and assay buffers do not contain free phosphate. This can be conveniently done by adding 200 µL of the Reagent to 40 µL sample solution. The blank OD values at 620 nm should be lower than 0.3. If the OD readings are higher than 0.3, check phosphate level. Lab detergents may contain high levels of phosphate. Make sure that lab wares are free from contaminating phosphate after thorough washes.

### ACTIVITY DETERMINATION IN 96-WELL PLATE

1. *Preparation of phosphate standards.* Prepare 800 µL Premix solution containing 50 µM phosphate by mixing 40 µL 1 mM phosphate standard with 760 µL distilled water. Number the tubes. Dilute standards as shown in the following Table. Pipette 40 µL standard in duplicate into wells of a clear-bottom 96-well plate.

No	Premix + H <sub>2</sub> O	Final Vol (µL)	Phosphate Conc (µM)	pmoles Phosphate in 40 µL
1	200µL + 0µL	200	50	2,000
2	160µL + 40µL	200	40	1,600
3	120µL + 80µL	200	30	1,200
4	80µL + 120µL	200	20	800
5	60µL + 140µL	200	15	600
6	40µL + 160µL	200	10	400
7	20µL + 180µL	200	5	200
8	0µL + 200µL	200	0	0

2. *Perform a series dilution* of enzyme in assay buffer. Set up 40-µL reactions and a control with no enzyme in separate wells. Incubate the reaction for desired period of time (e.g. 30 min).

Reaction Well

20 µL Assay Buffer  
10 µL enzyme  
10 µL 4 mM ATP or GTP

Control Well

20 µL Assay Buffer  
10 µL H<sub>2</sub>O  
10 µL 4 mM ATP or GTP

#### Gentaur Europe

Avenue de l'Armée 68  
B-1040 Brussels  
Tel: +32 16 58 90 45  
Fax: +32 16 50 90 45  
info@gentaur.com  
www.gentaur.com

#### Gentaur Germany

IME, Forckenbeckstrasse 6  
D-52074 Aachen Germany  
Tel: 0241 6085 13140  
Fax: 0241 6085 33033  
info@bioprice.com  
www.bioprice.com

#### Gentaur France

9, Rue Lagrange  
75005 Paris, France  
Tel: +33(0)143 25 0150  
Fax: +33(0)143 25 0160  
info@clonagen.com  
www.clonagen.com

#### Gentaur Italy

Tel: 02 36 00 65 93  
Fax: 02 36 00 65 94  
info@labprice.com  
www.labprice.com



[www.clonagen.com](http://www.clonagen.com)

[www.bioxys.com](http://www.bioxys.com)

[www.gentaur.com](http://www.gentaur.com)

3. Add 200  $\mu$ L Reagent and incubate 30 min at room temperature. Please note: use of a multi-channel pipettor is recommended. The Reagent terminates the enzyme reaction and generates color with the free phosphate produced in the enzyme reaction.

4. Read OD<sub>620nm</sub> on a plate reader.

5. *Enzyme activity.* Calculate  $\Delta$ OD values by subtracting OD values in reaction and control wells. Choose an enzyme concentration that gives a  $\Delta$ OD of 0.5 to 1, this will ensure that substrate hydrolysis (<10%) is within the linear kinetics of reaction. Compute the concentration of free phosphate produced from the standard curve. 1 Unit of activity is the amount of enzyme that catalyzes the production of 1  $\mu$  mole of free phosphate per minute under the assay conditions.

### INHIBITOR ASSAY IN 96-WELL PLATE

To evaluate an inhibitor or perform HTS, use the optimal enzyme concentration determined above. Incubate enzyme and inhibitor first for a certain period of time, before adding the substrate. At the end of reaction, add 200  $\mu$ L Reagent for phosphate determination.

#### Reaction Well

20  $\mu$ L Assay Buffer

5  $\mu$ L enzyme

5  $\mu$ L inhibitor

10  $\mu$ L 4 mM ATP or GTP

#### Control Well

20  $\mu$ L Assay Buffer

10  $\mu$ L 4 mM ATP or GTP

10  $\mu$ L Buffer/DMSO

### ASSAYS IN 384-WELL PLATE

The procedure is similar as in the 96-well plate assay, except that 20  $\mu$ L standards or 20  $\mu$ L reaction mixture (10  $\mu$ L Assay Buffer, 5  $\mu$ L 4 mM ATP, 5  $\mu$ L enzyme) are mixed with 80  $\mu$ L Reagent.

### GENERAL CONSIDERATIONS

*Materials.* Use ultrapure ATP and GTP. The provided 2x assay buffer contains 40 mM Tris, 80 mM NaCl, 8 mM MgAc<sub>2</sub>, 1 mM EDTA, pH 7.5. Other buffers (Hepes, Mes, Mops) can be used. Assay is compatible with 1 mM DTT, 2mM  $\beta$ -mercaptoethanol, 0.5 mg/mL BSA and 5% DMSO.

#### Gentaur Europe

Avenue de l'Armée 68  
B-1040 Brussels  
Tel: +32 16 58 90 45  
Fax: +32 16 50 90 45  
info@gentaur.com  
www.gentaur.com

#### Gentaur Germany

IME, Forckenbeckstrasse 6  
D-52074 Aachen Germany  
Tel: 0241 6085 13140  
Fax: 0241 6085 33033  
info@bioprice.com  
www.bioprice.com

#### Gentaur France

9, Rue Lagrange  
75005 Paris, France  
Tel: +33(0)143 25 0150  
Fax: +33(0)143 25 0160  
info@clonagen.com  
www.clonagen.com

#### Gentaur Italy

Tel: 02 36 00 65 93  
Fax: 02 36 00 65 94  
info@labprice.com  
www.labprice.com