

FavorPrep™

96-Well Viral DNA/ RNA Extraction Kit

Cat. No. FAVRE 001 (4 Plates)

FAVRE 002 (10 Plates)

Store at room temperature (15 - 25°C)

For research use only

Introduction

FavorPrep 96-Well Viral DNA/ RNA Kit is designed for high-throughput extraction of viral DNA/ RNA from cell free sample such as serum, plasma, body fluids and the supernatant of cell cultures. The method use a specialized chaotropic salt to lyse viruses, then DNA/ RNA in chaotropic salt is bonded to glass fiber matrix when the lysis mixture passing through the binding plate. After washing off the contaminants, the purified DNA/ RNA is eluted by RNase free water. The entire procedure can be completed in one hour without phenol/ chloroform extraction and alcohol precipitation. In the procedure, DNA/ RNA binding and washing steps could be do on vacuum manifold or by centrifuge. And the eluted DNA RNA is ready to use in different application.

Quality Control

The quality of 96-Well Viral DNA/ RNA Kit is tested on a lot-to-lot basis. The purified RNA is checked by agarose gel analysis and quantified with spectrophotometer.

Caution

VNE Buffer and Wash 1 Buffer contain irritant agent. Wear gloves and lab coat when handling these buffer.

Kit Contents

Name	FATRE001 4 plates	FATRE002 10 plates
VNE Buffer	200 ml	500 ml
AD Buffer (concentrated)*	26 ml	26ml X2
Wash 1 Buffer	130 ml	130 ml x 3
Wash 2 Buffer (concentrated)*	50 ml	50 ml x 3
RNase-free ddH ₂ O	30ml	60 ml
96-Well viral RNA/ DNA binding plate	4 pcs	10 pcs
96-Well 0.35 µl collection plate	4 pcs	10 pcs
Adhesive film	8 pcs	20 pcs

*Add 200 ml of ethanol (96-100%) to AD Buffer and Wash 2 Buffer when first use.

References

- (1) Vogelstein, B., and Gillespie, D. (1979) Proc. Natl. Acad. Sci. USA 76, 615.

Centrifuge Protocol

Step 1

lysis

- Transfer **200 µl of sample** to each well of the 96-Well 2 ml plate. If prepared sample is less than 200 µl, adjust sample volume to 200 µl with PBS (not provided)
- Add **400 µl of VNE Buffer** to each well of the 96-Well 2 ml plate.
- Seal the plate with Adhesive Film and mix by vortex. And Incubate at room temperature for 10 minutes.

Step2

Binding

- Add **450 µl of AD Buffer (ethanol added)** to each well of the 96-Well 2 ml plate. Mix the sample completely by shaking (seal with adhesive film) or pipetting.
- Place a **96-Well Viral DNA/ RNA Binding Plate** on top of another 96-Well 2 ml plate (not provided).
- Transfer the sample mixture to each well of the **96-Well Viral DNA/ RNA Binding Plate**.
- Place the combined plates in a rotor bucket and centrifuge at 3,500 rpm (1,000 ~1,500 x g) for 5 minutes.
- Discard the flow-through and place the **96-Well RNA Binding Plate** back to top of the 96-Well 2 ml plate.

Step 3

Wash

- Add **300 µl of Wash 1 Buffer** to each well of the **96-Well RNA Binding Plate**.
- Place the combined plates in a rotor bucket and centrifuge at 3,500 rpm (1,000 ~1,500 x g) for 5 minutes.
- Add **600 µl of Wash 2 Buffer (ethanol added)** to each well of the **96-Well Viral DNA/ RNA Binding Plate**.
- Place the combined plates in a rotor bucket and centrifuge at 3,500 rpm (1,000 ~1,500 x g) for 5 minutes.
- Discard the flow-through and place the **96-Well Viral DNA/ RNA Binding Plate** back to top of the 96-Well 2 ml plate.
- Place the combined plates in a rotor bucket and centrifuge at 3,500 rpm (1,000 ~1,500 x g) for an additional 10 minutes to remove residual ethanol.

Step 4

Elution

- Transfer the **96-Well Viral DNA/ RNA Binding Plate** on a clean 96-Well 350µl collection plate. (provided)
- Add **50 µl** RNase-free ddH₂O in the center of each well of DNA/ RNA Binding Plate.
- Stand for 3 minutes until RNase-free ddH₂O absorbed by the matrix. Centrifuge for 5 min at 3,500 rpm to elute purified RNA into the 96-Well 350µl collection plate. Seal with Adhesive Film and store 4°C or -20°C

Vacuum/ Centrifuge Protocol

Step 1

lysis

- Transfer **200 µl of sample** to each well of the 96-Well 2 ml plate. If prepared sample is less than 200 µl, adjust sample volume to 200 µl with PBS (not provided)
- Add **400 µl of VNE Buffer** to each well of the 96-Well 2 ml plate.
- Seal the plate with Adhesive Film and mix by vortex. And Incubate at room temperature for 10 minutes.

Step 2

Binding

- Add **450 µl of AD Buffer (ethanol added)** to each well of the 96-Well 2 ml plate. Mix the sample completely by shaking (seal with adhesive film) or pipetting.
- Place a **96-Well Viral DNA/ RNA Binding Plate** on top of the vacuum manifold.
(optional) Place a 96-Well 2 ml plate inside to collect waste.
- Transfer the sample mixture to each well of the **96-Well Viral DNA/ RNA Binding Plate**.
- Apply vacuum at 10 inches Hg for 5 minutes until wells have emptied.

Step 3

Wash

- Add **300 µl of Wash 1 Buffer** to each well of the **96-Well Viral DNA/ RNA Binding Plate**.
- Apply vacuum at 10 inches Hg for 5 minutes until wells have emptied.
- Add **600 µl of Wash 2 Buffer (ethanol added)** to each well of the **96-Well Viral DNA/ RNA Binding Plate**.
- Apply vacuum at 10 inches Hg for 5 minutes until wells have emptied.
- Apply vacuum at 10 inches Hg for an additional 10 minutes (or incubate at 60 °C for 10 minutes) to remove residual ethanol.

Step 4

Elution

- Transfer the **96-Well Viral DNA/ RNA Binding Plate** on a clean 96-Well 350µl collection plate. (provided)
- Add **50 µl of RNase-free ddH₂O** in the center of each well of RNA Binding Plate.
- Stand for 3 minutes until RNase-free ddH₂O absorbed by the matrix.
- Place the combined plates in a rotor bucket and centrifuge at 3,500 rpm for 5 min to elute purified DNA/ RNA into the 96-Well 350µl collection plate. Seal with Adhesive Film and store 4°C or -20°C.