

ELISA Kit for IgM Antibody to Hepatitis A Virus

Catalog No.: KO31008096

[NAME AND INTENDED USE]

ELISA Kit for IgM Antibody to Hepatitis A Virus is an *in vitro* enzyme immunoassay for the detection of HAV-IgM in human serum or plasma.

[PRINCIPLE]

This kit uses capture ELISA method to detect anti-HAV IgM. The purified rabbit anti human IgM monoclonal antibody (Anti- μ -chain) is coated on the solid phase of multi-wells. Serum sample, Horseradish peroxidase labeled HAVAg are added to coated wells. After incubation, if HAV-IgM is present in the sample, a complex of Anti- μ -chain-HAV-IgM-HAVAg labeled with HRP will form. Wash wells to remove other unbound serum components, incubate with substrate (TMB) to form a colored product, and measure the absorbance at 450 nm to indicate the presence or absence of HAV-IgM in the sample. The test is special, sensitive, reproducible and easy to operate.

[MATERIALS PROVIDED]

1. Anti- μ -chain Coated Microwell Plate	1 block (96wells)
2. Enzyme Conjugant (HAVAg-HRP)	1 bottle (12ml)
3. Negative Control Serum	1 vial (1ml)
4. Positive Control Serum	1 vial (1ml)
5. Wash Buffer (1:20 dilution prior to use)	1 bottle (50ml)
6. Substrate A	1 bottle (6ml)
7. Substrate B	1 bottle (6ml)
8. Stop Solution	1 bottle (6ml)
9. Plastic Bag	1 bag
10. Seal Paper	3 pieces
11. Manual	1 each

[SAMPLE COLLECTION AND PRESERVATION]

Blood serum samples are routinely prepared from vein. Blood plasma sample are routinely prepared with routine amount of anticoagulant such as heparin or sodium citrate. Sample can be stored at 4°C if tested within five days. Sample can be stored at -20°C at least for 3 months. Avoid hemolysis and repetitive freeze and thaw of samples. Samples with cloud or precipitation should be centrifugated or filtered before test. Prevent serum from bacteria contamination during collection and storage.

[TEST PROCEDURE]

- Bring *** ELISA Kit for IgM Antibody to Hepatitis A Virus *** (all reagents), and samples to room temperature before use (approximately 30 minutes).
- Dilute concentrated wash buffer 1:20 with ddH₂O
- For each test, set one blank, two positive and three negative controls. Add 50 μ l positive and negative control serum into positive and negative control wells respectively without sample diluent.
- Add 50 μ l test serum into test wells.
- Cover wells with seal paper, incubate for 30 minutes at 37°C.
- Discard the liquid in all wells and fill the wells with wash solution (300 μ l per well). Lay aside for 15 seconds, discard the liquid in all wells and fill the wells with wash solution. Repeat 5 times and dry wells after last wash.
- Add 100 μ l Enzyme Conjugant in each well except the blank well.
- Cover wells with seal paper, incubate for 30 minutes at 37°C.
- Wash 5 times the same as step 6.
- Add 50 μ l substrate A and B respectively to each well,

- mix gently, protected from light and lay aside for 15 minutes at 37°C.
- Add one drop of stop solution (50 μ l) into each well to stop the reaction, including blank well.
 - Measure the absorbance at 450 nm against the blank, or measure the absorbance at 450 nm/630-690 nm.

[INTERPRETATION OF RESULTS]

Colorimetric Method

Cut Off Value calculation:

COV = 0.07 + the average OD of negative controls

Positive OD₄₅₀ of sample \geq COV

Negative OD₄₅₀ of sample < COV

Notes If the absorbance of negative controls is below 0.05, calculate it as 0.05. If the absorbance of negative controls is above 0.05, calculate it as its original value.

[PRECAUTIONS]

- The samples should be fresh, avoid hemolysis, bacteria growing, and repetitive freeze and thaw.
- Do not interchange reagents between kit lots. The seal paper can't be used repeatedly.
- Mix reagents well before use. If crystal form in certain reagents, such as wash buffer, it can be used without problems after warm up and mix well.
- Follow instruction exactly during assay, especially in temperature and time for reactions. All pipetting devices should be used with care and calibrated regularly following the manufacturer's instructions.
- Put the remained reagents to the sealed pouch, and return to 2~8°C in time.
- To prevent cross-contamination, wear gloves and working suits throughout the procedure, and execute the disinfection and isolation regulations strictly. Dispose of all samples and materials used to perform the test. The 5.0g/L liquid sodium hypochlorite solution or 121°C high pressure steam may be used to disinfect samples and materials before disposal (The positive control serum in the kit has been inactivated already)

[PACKAGE SIZE]

96 tests/Kit

[STORAGE AND STABILITY]

Store the kit at 2~8°C.

[EXPIRATION]

The shelf life is 12 months from the receiving date.

This Kit is for Research Use Only