


Expression Arrest™ GIPZ lentiviral shRNAmir library





openbiosystems.com info@openbiosystems.com
6705 Odyssey Drive /P/ 888.412.2225
Huntsville, AL 35806 /F/ 256.704.4849

Expression Arrest™ Human GIPZ lentiviral shRNAmir library RHS4430, RHS4477

The human GIPZ lentiviral shRNAmir library was developed by Open Biosystems in collaboration with Dr. Greg Hannon (CSHL) and Dr. Steve Elledge (Harvard). This library combines the design advantages of microRNA-adapted shRNA (shRNAmir) with the pGIPZ lentiviral vector to create a powerful RNAi trigger capable of producing RNAi in most cell types including primary and non-dividing cells.

shRNAmir triggers have been designed to mimic a natural microRNA primary transcript and each target sequence has been selected based on thermodynamic criteria for optimal small RNA performance. Validation of this design is detailed in Silva *et al* (2005) showing a substantial increase in knockdown efficiency.

Unique features of the GIPZ lentiviral shRNAmir library include:

- shRNAmir constructs targeting the entire human genome already cloned into the pGIPZ lentiviral vector
- Efficient low copy knockdown - Important for pooled screens
- TurboGFP (tGFP) and shRNAmir are part of a bicistronic transcript allowing the visual marking of shRNAmir expressing cells
- Effective transduction of primary and non-dividing cell lines e.g. neurons
- Unique 60nt molecular barcode facilitate pooled screens

Shipping and Storage

The Expression Arrest Human GIPZ shRNAmir lentiviral library is provided in 96-well microtiter plates containing frozen stock cultures of *E. coli* (Prime+) in LB-Lennox (low salt) broth with 8% glycerol, 100µg/ml carbenicillin and 25µg/ml zeocin.

Individual constructs are shipped as bacterial cultures of *E. coli* (prime+) in LB-Lennox (low salt) broth with 8% glycerol and carbenicillin (100ug/ml) and zeocin (25ug/ml). Individual constructs are shipped on wet ice. Open Biosystems checks all cultures for growth prior to shipment.

The GIPZ human lentiviral shRNAmir library and individual constructs should be stored at -80°C.




Figure 2: pGIPZ lentiviral vector

Table 1: Features of the pGIPZ Vector

Vector Element	Utility
CMV Promoter	RNA Polymerase II promoter
cPPT	Central Polypurine tract helps translocation into the nucleus of non-dividing cells
WRE	Enhances the stability and translation of transcripts
tGFP	Marker to track shRNAmir expression
IRES-Puro	Mammalian selectable marker
AMP ^r	Ampicillin bacterial selectable marker.
5'LTR	5' long terminal repeat
pUC ori	High copy replication and maintenance of plasmid in <i>E.coli</i>
SIN-LTR	3' Self inactivating long terminal repeat
RRE	Rev response element
ZEO ^r	Bacterial selectable marker

Table 2: Antibiotic Resistances Conveyed by pGIPZ

Antibiotic	Concentration	Utility
Ampicillin (carbenicillin)	100µg/ml	Bacterial selection marker
Zeocin	25µg/ml	Bacterial selection marker (vector)
Puromycin	variable	Mammalian selectable marker

Culturing protocols and maintenance of pGIPZ

It is well known that viral vectors have a tendency to recombine producing background recombinants. Recombination occurs at the long terminal repeat regions (LTR's). The LTR recombination, which results in loss of most of the plasmid, can confer a growth advantage on the cells. It is therefore critical to maintain careful growth conditions when culturing viral vectors in *E.coli* in order to reduce the number and abundance of background recombinants. The GIPZ lentiviral shRNAmir library has passed through internal QC processes to ensure high quality and low recombination.




Figure 3. Representative shRNAmir containing pGIPZ lentiviral clones grown for 16 hours at 30° C and the plasmid isolated and normalized to a standard concentration. Clones were then digested with *Sac*II and run out on a gel. The expected band sizes are (bp)= **1259, 2502, 7927**. No recombinant products are visible. 10kb molecular weight ladder (10kb, 7kb, 5kb, 4kb, 3kb, 2.5kb, 2kb, 1.5kb, 1kb)




Figure 4. Gel image of a single plate from the GIPZ library cultured for 10 successive generations in an attempt to determine the tendency of the pGIPZ vector to recombine. Each generation was thawed, replicated and incubated O/N for 16 hours at 30° C then frozen, thawed and replicated. This process was repeated for 10 growth cycles. After the 10th growth cycle, plasmid was isolated and normalized to a standard concentration. Clones were then digested with *Sac*II and run on a gel. Expected band sizes (bp) = **1259, 2502, 7927**. 10kb molecular weight ladder (10kb, 7kb, 5kb, 4kb, 3kb, 2.5kb, 2kb, 1.5kb, 1kb) The pGIPZ vector appears stable without showing any recombination.

Culture conditions for individual plasmid preparations

Most plasmid mini-prep kits recommend a culture volume of 1–10ml for good yield. For shRNAmir constructs, 5ml of culture can be used for one plasmid mini-prep generally producing 5–10µg of plasmid DNA.

1. Upon receiving your glycerol stock(s) containing the shRNAmir of interest store at –80°C until ready to begin.
2. To prepare plasmid DNA first thaw your glycerol stock culture and pulse vortex to resuspend any *E. coli* that may have settled to the bottom of the tube.
3. Take a 10µl inoculum from the glycerol stock into 3-5ml of LB (low salt) with 100µg/ml carbenicillin and 25ug/ml zeocin. Incubate at 37°C for 16 hours with vigorous shaking. Return the glycerol stock(s) to -80°C. If a larger culture volume is desired, use the 3-5ml overnight culture as a starter inoculum. Incubate at 37°C for 16 hrs with vigorous shaking.
4. Pellet the 3-5ml culture and begin preparation of plasmid DNA.
5. Run 3-5µl of the plasmid DNA on a 1% agarose gel. pGIPZ with shRNAmir is 11744bp.

Note: Due to the tendency of all viral vectors to recombine we recommend keeping the incubation times as short as possible and avoid subculturing. Return to your original glycerol stock or the colony glycerol stock for each plasmid preparation.

Restriction Digests of pGIPZ

The following is a sample protocol for restriction enzyme digestion using *KpnI*, *SacII*, *Sall*, *XhoI* and/or *NotI* for diagnostic quality control of pGIPZ lentiviral vectors.

1. Using filtered pipette tips and sterile conditions add the following components, in the order stated, to a sterile PCR thin-wall tube.

Sterile, nuclease-free water	X µl
Restriction enzyme 10X buffer	1µl
BSA (10X, 10mg/ml) if required	1µl
DNA sample 80 -240ng, in water or TE buffer	X µl
<u>Restriction enzyme 20U</u>	<u>0.25µl</u>
Final volume	10µl

2. Mix gently by pipetting.
3. Incubate in a thermalcycler at 37°C for 2 hours to digest
4. Load the gel with 10µl of each of the digested samples (*KpnI*, *SacII*, *Sall*, *XhoI* and/or *NotI*) on a 1% agarose gel. Run uncut sample alongside the digested samples.




Figure 5: Restriction digests with pGIPZ. Lane 1– 10kb molecular weight ladder (10kb, 7kb, 5kb, 4kb, 3kb, 2.5kb, 2kb, 1.5kb, 1kb). Lane2 - Uncut pGIPZ vector. Lane 3 - *KpnI* digested pGIPZ produces 2 bands at 1750bp and 9860bp. Lane4- *SacII* digest produces 3 bands at 1178bp, 2502bp and 7930bp. Lane 5 -*SalI* produces 3 bands at 2188bp, 4298bp and 5124bp. Lane 6 – *XhoI NotI* double digest produces 2 bands at 1210bp and 10400bp.

Culture conditions for 96-well plasmid preparation

Inoculate 96-well bio-block containing 1ml per well of the above media with 1µl of the culture. Incubate at 37°C with shaking (~170-200 RPM). We have observed that incubation times from 16 hours produces good plasmid yield. For plasmid preparation, follow the kit protocols recommended by the manufacturer.

Note: The cells can be grown at 37°C for purposes of template preparation or sequencing. For archive replication, grow all pGIPZ clones at 30°C in LB-Lennox (low salt) media plus 25ug/ml zeocin and 100ug/ml carbenicillin in order to provide maximum stability of the clones.

Materials Required

LB-Lennox Broth (low salt) – VWR item# EM1.00547.0500
 Glycerol – VWR item# EM-4760
 Carbenicillin or Ampicillin – VWR item# EM-2200 or 80030-956
 Zeocin – Invivogen item# ant-zn-5p
 96-well microplates – VWR item# 62407-174
 Aluminum seals – VWR item# 73520-056
 Disposable replicators – Genetix item# X5054

CaPO₄ Transfection Protocol for pGIPZ Lentiviral packaging

(100-mm dish format)

1. Approximately 24 hours before transfection, seed 6.0×10^6 293T cells in 14-ml of complete media (Dulbecco's modified Eagle's medium (DMEM), 10% FBS, 2 mM L-glutamine, 1X Pen-Strep).
2. Incubate at 37°C, 5% CO₂ overnight. Transfection should begin when cells are approximately 90% confluent.
3. The following describes the preparation of DNA-CaPO₄ mixture and the protocol for performing 1 transfection (one 100-mm dish). Transfection reactions are carried out in a 5-ml polystyrene round-bottom tube (Falcon catalog # 352058)

DNA Preparation

DNA to be co-transfected, add volume to 945 µl with sterile water:

- | | |
|--|--------|
| 1. Transgene (gene transfer vector): | 21µg |
| 2. pCMV-Gag-Pol (2 nd Generation eg. psPAX2): | 21µg |
| 3. pCMV-VSV-G-poly A (e.g pMD2.G) | 10.5µg |

Note: The number of transfection reactions is scalable. For example, if transfecting numerous 100-mm dishes to generate larger volumes or higher titers of the same vector stock, a master mix of the DNA-water stock is made and aliquoted into 50-ml **polystyrene** tubes. A maximum of seven 100-mm dishes can be transfected from one 50-ml tube. For seven transfections, pipette 6615 µl of the DNA-water mix maintaining the same ratio of each of the vector plasmids as well as DNA to water.

4. The following describes the CaPO₄ precipitation reaction in both one and seven 100-mm dish formats.

One 100-mm dish:

In one 5-ml snap cap polystyrene tube mix:

- a. DNA plus sterile water to final volume of 945 µl.
- b. Add 105 µl of 2.5 M CaCl₂.
- c. While vortexing tube, add dropwise 1050 µl of 2X HBSS (2100 µl total volume). Make sure vortexer is set so that the contents mix thoroughly without spilling over.

For seven 100-mm dishes:

In one 50-ml polystyrene tube mix:

- a. DNA plus sterile water to final volume of 6615 µl.
- b. Add 735 µl of 2.5 M CaCl₂.
- c. While vortexing tube, add dropwise 7350 µl of 2X HBSS (14,700 µl total volume). Make sure vortexer is set so that the contents mix thoroughly without spilling over.

5. Incubate at room temperature for 3 minutes. A chalky white precipitate should be visible in the tube. If no precipitate is noticeable, allow the incubation to continue at room temperature until it is visible.
6. Following incubation, vortex contents of the tube a few seconds, and pipette 2100 μ l of the transfection mixture dropwise into one well. Do not add the transfection mixture to only one area of the well but instead spread the drops over the entire surface of well.
7. Incubate at 37°C, 5% CO₂ for 12-16 hours.
8. Remove media from each plate and slowly pipette 14 ml of DMEM, **5% FBS**, 2 mM L-glutamine, 1X Pen-Strep) to each well. **DO NOT WASH** cells. 5% FBS is used to decrease the amount of serum proteins pelleted with the Vector stock during ultracentrifugation.
9. Incubate at 37°C, 5% CO₂ for an additional 48 hours.
10. Harvest virus-containing supernatant. Pellet cells/debris by low-speed centrifugation (1600 x g for 10 min).
11. Aliquot virus and store at -80°C.
12. Virus can be concentrated by ultracentrifugation (SW28, 23,000rpm, 1.5h @ 4°C).

Reagents:

2.5 M CaCl₂

(For 100 ml):

36.75 g CaCl₂ (Sigma, Cat. No. C-7902)
Add sterile dH₂O to 100 ml

Filter-sterilize through 0.22 μ m filter flask (Millipore)

2X HBSS (Hepes Buffered Saline Solution)

50 mM Hepes (pH 7.1)
280 mM NaCl
1.5 mM Sodium Phosphate

The final pH should be 7.1

(For 1 liter):

11.915 g	Hepes (Sigma, Cat. No. H-3375)
16.363 g	NaCl (Sigma, Cat. No. S-3014)
0.090 g	NaH ₂ PO ₄ (Sigma, Cat. No. S-3139)
0.107 g	Na ₂ HPO ₄ (Sigma, Cat. No. S-3264)

Add sterile dH₂O to 990 ml
pH to 7.1 by dropwise adding 10 N NaOH




Figure 6: Detailed Vector Map of pGIPZ lentiviral vector

Sequence of pGIPZ lentiviral vector (11774bp)

```

5'LTR(Lenti-WT) other(1,635)>>>
|
U3(HIV-LTR) reg(1,455)>>>
|
1  tgggaagggctaattcactcccaagaagacaagatatccttgatctgtggatctaccaca 60
   ACCTTCCCATTAAAGTGAGGGTTTCTTCTGTTCTATAGGAAGTACACCTAGATGGTGT

61  cacaaggctacttccttgattagcagaactacacaccagggccaggggtcagatatccac 120
   GTGTTCCGATGAAGGGACTAATCGTCTTGATGTGTGGTCCCGGTCCCCAGTCTATAGGTG

121  tgacctttggatggtgctacaagctagtaccagttgagccagataaggtagaagaggcca 180
   ACTGGAAACCTACCACGATGTTTCGATCATGGTCAACTCGGTCTATTCCATCTTCTCCGGT

181  ataaaggagagaacaccagcttgttacaccctgtgagcctgcatgggatggatgaccggy 240
   TATTTCTCTCTTGTGGTTCGAACAATGTGGGACACTCGGACGTACCCTACCTACTGGGCC

```

241 agagagaagtgttagagtggaggtttgacagccgcctagcatttcatcacgtggcccgag 300
TCTCTCTTCACAATCTCACCTCCAAACTGTGCGCGGATCGTAAAGTAGTGCACCGGGCTC

301 agctgcatccggagtacttcaagaactgctgatatcgagcttgctacaagggactttccg 360
TCGACGTAGGCCTCATGAAGTTCTTGACGACTATAGCTCGAACGATGTTCCCTGAAAGGC

361 ctggggactttccagggagggcgtggcctgggaggactggggagtgggcagccctcagat 420
GACCCCTGAAAGGTCCCTCCGCACCGGACCCGCCCTGACCCCTCACCGCTCGGGAGTCTA

R(HIV-LTR) reg(456,550)>>>
|

421 cctgcatataagcagctgctttttgacctgactgggtctctctggttagaccagatctga 480
GGACGTATATTCGTGACGAAAAACGGACATGACCCAGAGAGACCAATCTGGTCTAGACT

481 gcctgggagctctctggctaactagggaaaccactgcttaagcctcaataaagcttgct 540
CGGACCCTCGAGAGACCGATTGATCCCTTGGGTGACGAATTCGGAGTTATTTTCGAACGGA

U5(HIV-LTR) reg(551,635)>>>
|

541 tgagtgttcaagtagtgtgtgcccgtctggtgtgactctggtaactagagatccctc 600
ACTCACGAAGTTCATCACACACGGGCAGACAACACACTGAGACCATTGATCTCTAGGGAG

601 agacccttttagtcagtggtgaaaatctctagcagtgccgcccgaacagggacttgaaag 660
TCTGGGAAAATCAGTCACACCTTTTAGAGATCGTCACCGCGGGCTTGTCCCTGAACTTTC

PSI(HIV) reg(685,822)>>>
|

661 cgaaagggaaaccagaggagctctctcgacgcaggactcggcttgctgaagcgcgcacgg 720
GCTTTCCTTTGGTCTCCTCGAGAGAGCTGCGTCCTGAGCCGAACGACTTCGCGCGTGCC

721 caagaggcgagggggcggcactgggtgagtacgccaaaaatttgactagcggaggctaga 780
GTTCTCCGCTCCCCGCGCTGACCACTCATGCGGTTTTTAAACTGATCGCCTCCGATCT

NruI
|

781 aggagagagatgggtgagagagcgtcagtattaagcgggggagaattagatcgcgatggg 840
TCCTCTCTTACCCAGCTCTCGCAGTCATAATTCGCCCCCTCTTAATCTAGCGCTACCC

841 aaaaaattcggttaaggccagggggaaagaaaaaataaaattaaaacatatagtatggg 900
TTTTTTAAGCCAATTCCGGTCCCCCTTTCTTTTTTATATTTAATTTTGTATATCATACCC

901 caagcagggagctagaacgattcgcagttaatcctggcctgtagaaacatcagaaggct 960
GTTTCGTCCTCGATCTTGCTAAGCGTCAATTAGGACCGACAATCTTTGTAGTCTTCCGA

961 gtagacaaatactgggacagctacaaccatcccttcagacaggatcagaagaacttagat 1020
CATCTGTTTTATGACCCTGTCGATGTTGGTAGGGAAGTCTGTCCTAGTCTTCTTGAATCTA

1021 cattatataatacagtagcaaccctctattgtgtgcatcaaaggatagagataaaagaca 1080
GTAATATATTATGTCATCGTTGGGAGATAACACACGTAGTTTCCTATCTCTATTTTCTGT

1081 ccaaggaagcttttagacaagatagaggaagagcaaaaacaaaagtaagaccaccgcacagc 1140
GGTTCCTTCGAAATCTGTTCTATCTCCTTCTCGTTTTTGTTCATTCTGGTGGCGTGTGCG

1141 aagcggccggccgctgatcttcagacctggaggaggagatatgagggacaattggagaag 1200
TTCGCCGGCCGGCGACTAGAAGTCTGGACCTCCTCCTCTATACTCCCTGTTAACCTCTTC

1201 tgaattatataaaatataaaagtagtaaaaaattgaaccattaggagtagcaccaccaaggc 1260
ACTTAATATATTTTATATTTTCATCATTTTTAACTTGGTAATCCTCATCGTGGGTGGTTCCG

RRE (HIV)

reg(1314,1518)>>>

1261 aaagagaagagtggtgcagagagaaaaaagagcagtggggaataggagctttgttccttgg 1320
TTTCTCTTCTCACCACGTCTCTTTTTTCTCGTCACCCTTATCCTCGAAACAAGGAACC

1321 gttcttgggagcagcaggaagcactatgggagcagcgtcaatgacgctgacggtagcagc 1380
CAAGAACCCTCGTCGTCTTTCGTGATACCCGCGTCGCAGTTACTGCGACTGCCATGTCCG

1381 cagacaattattgtctggtatagtgagcagcagacaatttgctgagggctattgagggc 1440
GTCTGTTAATAACAGACCATATCACGTGTCGTCTTGTAAACGACTCCCATAACTCCG

1441 gcaacagcatctggtgcaactcacagctctggggcatcaagcagctccaggcaagaatcct 1500
CGTTGTCGTAGACAACGTTGAGTGTGACACCCCGTAGTTCGTGAGGTCCGTTCTTAGGA

1501 ggctgtggaaagatacctaaaggatcaacagctcctggggatttgggggtgctctggaaa 1560
CCGACACCTTTCTATGGATTTCTAGTTGTGAGGACCCCTAAACCCCAACGAGACCTTT

1561 actcatttgcaccactgctgtgccttgggaatgctagttggagtaataaatctctggaaca 1620
TGAGTAAACGTGGTGACGACACGGAACCTTACGATCAACCTCATTATTTAGAGACCTTGT

1621 gatttgggaatcacacgacctggatggagtgggacagagaaattaacaattacacaagctt 1680
CTAAACCTTAGTGTGCTGGACCTACCTCACCTGTCTCTTTAATTGTTAATGTGTTTCGAA

1681 aatacactccttaattgaagaatcgcaaaaccagcaagaaaagaatgaacaagaattatt 1740
TTATGTGAGGAATTAACCTTCTTAGCGTTTTGGTCGTTCTTTTCTTACTTGTCTTAATAA

1741 ggaattagataaatgggcaagtttgggaattgggttaacataacaaattggctgtggta 1800
CCTTAATCTATTTACCCGTTCAAACACCTTAACCAAATTGTATTGTTTAAACCGACACCAT

1801 tataaaaattattcataatgatagtaggaggttggtaggtttaagaatagtttttgcgtgt 1860
ATATTTTAATAAGTATTACTATCATCCTCCGAACCATCAAATTCTTATCAAAAACGACA

1861 actttctatagtgaaatagagtttaggcagggatattcaccattatcgtttcagaccacct 1920
TGAAAGATATCACTTATCTCAATCCGTCCCTATAAGTGGTAATAGCAAAGTCTGGGTGGA

1921 cccaaccccgaggggacccgacaggcccgaaggaatagaagaagaaggtggagagagaga 1980
GGGTTGGGGCTCCCCTGGGCTGTCCGGGCTTCCTTATCTTCTTCTTCCACCTCTCTCTCT

1981 cagagacagatccattcgattagtgaaacggatcggcactgctgagcccaattctgcagac 2040
GTCTCTGTCTAGGTAAGCTAATCACTTGCCTAGCCGTGACGCACGCGGTTAAGACGTCTG

CTS reg(2064,2214)>>>

2041 aaatggcagatttcatccacaattttaaaagaaaaggggggattgggggggtacagtgag 2100
TTTACCGTCATAAGTAGGTGTTAAAATTTTCTTTTCCCCCTAACCCCATGTCACGTC

2101 gggaaagaatagtagacataatagcaacagacatacaaaactaaagaattacaaaaacaaa 2160
CCCTTTCTTATCATCTGTATTATCGTTGTCTGTATGTTTGATTTCTTAATGTTTTTGT

2161 ttacaaaaattcaaaattttcgggtttattacagggacagcagagatccagtttggttag 2220
AATGTTTTTAAGTTTTAAAGCCCAATAATGTCCCTGTGCTCTCTAGGTCAAACCAATC

ZeoR marker(2245,2619)<<<

2221 taccgggcccgcctctagtcggaatcagtcctgctcctcggccacgaagtgcacgcagtt 2280
ATGGCCCGGGCGAGATCAGGCCTTAGTCAGGACGAGGAGCCGGTGCTTCACGTGCGTCAA
2281 gccggccgggtcgcgccagggcgaaactcccgcacccacggctgctcgccgatctcggtcat 2340
CGGCCGGGCCAGCGCGTCCCCTTGAGGGCGGGGGTGCCGACGAGCGGCTAGAGCCAGTA
2341 ggccggcccgaggcgtcccggaagtctcgtggacacgacctccgaccactcggcgtaacag 2400
CCGGCCGGGCCTCCGCAGGGCCTTCAAGCACCTGTGCTGGAGGCTGGTGAGCCGCATGTC
2401 ctcgtccaggccgcgccacccacacccaggccaggggtggtgctcggcaccacctggctcctg 2460
GAGCAGGTCCGGCGCGTGGGTGTGGGTCCGGTCCCACAACAGGCCGTGGTGGACCAGGAC
2461 gaccgcgctgatgaacagggtcacgctcgtcccggaccacaccggcgaaagtgcctccac 2520
CTGGCGCGACTACTTGTCCCAGTGCAGCAGGGCCTGGTGTGGCCGCTTCAGCAGGAGGTG
2521 gaagtcccgggagaacccgagccgggtcgggtccagaactcgaccgctccggcgacgtcgcg 2580
CTTCAGGGCCCTCTTGGGCTCGGCCAGCCAGGTCTTGAGCTGGCGAGGCCGCTGCAGCGC

SfiI

EM7 prom(2620,2683)<<<

2581 cgcggtgagcaccggaacggcactgggtcaacttggccatgggtggccctcctatagtgagt 2640
GCGCCACTCGTGGCCTTGCCGTGACCAGTTGAACCGGTACCACCGGGAGGATATCACTCA
2641 cgtattatactatgccgatatactatgccgatgattaattgtcaacacgtgctgcaggtc 2700
GCATAATATGATACGGCTATATGATACGGCTACTAATTAACAGTTGTGCACGACGTCCAG

XbaI

CMV-IE-Promoter-Enhancer

prom(2738,3311)>>>
2701 cgaggttctagacgtattaccgccatgcattagttattaatagtaatacaattacggggtc 2760
GCTCCAAGATCTGCATAATGGCGGTACGTAATCAATAATTATCATTAGTTAATGCCCCAG
2761 attagttcatagcccatatatggagttccgcttacataacttacggtaaatggcccggc 2820
TAATCAAGTATCGGGTATATACCTCAAGGCGCAATGTATTGAATGCCATTTACCGGGCGG
2821 tggctgaccgcccacgacccccgcccattgacgtcaataatgacgtatggtcccatagt 2880
ACCGACTGGCGGGTTGCTGGGGGCGGGTAACTGCAGTTATTACTGCATACAAGGGTATCA
2881 aacgccaatagggactttccattgacgtcaatgggtggagatatttacggtaaaactgccc 2940
TTGCGGTTATCCCTGAAAGGTAACCTGCAGTTACCCACCTCATAAATGCCATTTGACGGGT
2941 cttggcagtacatcaagtgtatcatatgcccaagtagccccctattgacgtcaatgacgg 3000
GAACCGTCATGTAGTTACATAGTATACGGTTCATGCGGGGATAACTGCAGTTACTGCC
3001 taaatggcccgcctggcattatgccagtagacattatgggactttcctacttggca 3060
ATTTACCGGGCGGACCGTAATACGGGTTCATGTACTGGAATACCCTGAAAGGATGAACCGT
3061 gtacatctacgtattagtcacgctattaccatgggtgatgagggttttggcagtagatcaa 3120
CATGTAGATGCATAATCAGTAGCGATAATGGTACCCTACGCCAAAACCGTCATGTAGTT
3121 tgggcgtggatagcgggttgactcaggggatttccaagtctccacccattgacgtcaa 3180

```

ACCCGCACCTATCGCCAAACTGAGTGCCCTAAAGGTTTACAGAGGTGGGGTAACTGCAGTT
3181 tgggagtttgTTTTGGCACCAAAATCAACGGGACTTTCCAAAATGTCGTAACAACCTCCGC 3240
ACCCTCAAACAAAACCGTGGTTTTAGTTGCCCTGAAAGGTTTTACAGCATTGTTGAGGGC

3241 cccattgacgcaaatgggCGGTAGGCGTGTACGGTGGGAGGTCTATATAAGCAGAGCTCG 3300
GGGTAACCTGCGTTTTACCCGCCATCCGCACATGCCACCCTCCAGATATATTTCGTCTCGAGC

3301 tttagtgaaccgtcagatcgctggagacgccatccacgctgttttgacctccatagaag 3360
AAATCACTTGGCAGTCTAGCGGACCTCTGCGGTAGGTGCGACAAAACCTGGAGGTATCTTC

                                turboGFP tag(3390,4088)>>>
                                |
3361 acaccgactctactagaggatctgccaccatggagagcgacgagagcggcctgcccgccca 3420
TGTGGCTGAGATGATCTCCTAGACGGTGGTACCTCTCGCTGCTCTCGCCGGACGGGGCGGT

3421 tggagatcgagtgccgcatcaccggcaccctgaacggcgtggagtccgagctggtggggcg 3480
ACCTCTAGCTCACGGCGTAGTGGCCGTGGGACTTGCCGCACCTCAAGCTCGACCACCCGC

3481 gCGGAGAGGGGcCCCCGAGCAGGGCGCATGACCAACAAGATGAAGAGCACCAAAGGGCG 3540
CGCCTCTCCCGTGGGGGCTCGTCCCGGCGTACTGGTTGTTCTACTTCTCGTGGTTTTCCGC

3541 ccctgaccttcagcccctacctgctgagccacgtgatgggctacggcttctaccacttcg 3600
GGGACTGGAAGTCGGGGATGGACGACTCGGTGCACTACCCGATGCCGAAGATGGTGAAGC

3601 gcacctaccccagcgggctacgagaacccttctctgcacgccatcaacaacggcgggtaca 3660
CGTGGATGGGGTCGCCGATGCTCTTGGGGAAGGACGTGCGGTAGTTGTTGCCGCCGATGT

3661 ccaacaccgcatcgagaagtagcaggacggcggcgctgctgcacgtgagcttcagctacc 3720
GGTTGTGGGCGTAGCTCTTCATGCTCCTGCCGCCGCACGACGTGCACTCGAAGTCGATGG

3721 gctacgaggccggcggcgtgatcggcgacttcaaggtgatgggacccggcttccccgagg 3780
CGATGCTCCGGCCGGCGCACTAGCCGCTGAAGTTCCTACTACCCGTGGCCGAAGGGGCTCC

3781 acagcgtgatcttcaccgacaagatcatccgcagcaacgccaccgtggagcacctgcacc 3840
TGTCGCACTAGAAGTGGCTGTTCTAGTAGGCGTCTGTTGCGGTGGCACCTCGTGGACGTGG

3841 ccatgggCGATAACGATCTGGATGGCAGCTTcaccgCACCTTcagcctgcgCGACGGCG 3900
GGTACCCGCTATTGCTAGACCTACCGTGAAGTGGGCGTGAAGTCCGACGCGCTGCCGC

3901 gctactacagctccgtgggtggacagccacatgcacttcaagagcgccatccaccccagca 3960
CGATGATGTCGAGGCACCACCTGTGCGGTGTACGTGAAGTTCCTCGCGGTAGGTGGGGTCTGT

3961 tcctgcagaacgggggccccatgttcgccttccgcccgtggaggaggatcacagcaaca 4020
AGGACGTCTTGCCCCCGGGGTACAAGCGGAAGGCGGGCGCACCTCCTCCTAGTGTCTGTTGT

4021 ccgagctgggcatcgtggagtaccagcagccttcaagaccccgatgcagatgcccgtg 4080
GGCTCGACCCGTAGCACCTCATGGTCTGTCGGAAGTTCCTGGGGCCTACGTCTACGGCCAC

                                NotI          IRES reg(4114,4689)>>>
                                |              |
4081 aagaataatgtacaagtagcggcgcaaaattccgccccctctccctccccccccctaacg 4140
TTCTTATTACATGTTTCATCGCCGGCGTTTAAGGCGGGGAGAGGGAGGGGGGGGGATTGC

4141 ttactggccgaagccgcttgggaataaggccgggtgtgcgtttgtctatatgttatTTTTCCA 4200
AATGACCGGCTTCGGCGAACCTTATTCCGGCCACACGCAAACAGATATACAATAAAAGGT

```

4201 ccatattgccgtcttttggcaatgtgagggcccgaaacctggccctgtcttcttgacga 4260
GGTATAACGGCAGAAAACCGTTACTTCCCAGGCTTTGGACCGGGACAGAAGAACTGCT

4261 gcattcctaggggtctttcccctctcgccaaaggaatgcaaggtctggtgaatgtcgtga 4320
CGTAAGGATCCCCAGAAAGGGGAGAGCGGTTTCCTTACGTTCCAGACAACCTTACAGCACT

4321 aggaagcagttcctctggaagcttcttgaagacaaacaacgtctgtagcgcaccctttgca 4380
TCCTTCGTCAAGGAGACCTTCGAAGAACTTCTGTTTTGTTGCAGACATCGCTGGGAAACGT

4381 ggcagcggaaacccccacctggcgacaggtgcctctgcgccaaaagccacgtgtataag 4440
CCGTGCGCTTGGGGGTGGACCGCTGTCCACGGAGACGCCGGTTTTTCGGTGCACATATTC

4441 atacacctgcaaagggcgacacacccagtgccacggtgtgagttggatagttgtggaaa 4500
TATGTGGACGTTTTCCGCCGTGTTGGGGTACGGTGCAACACTCAACCTATCAACACCTTT

4501 gaggcaaatggctctcctcaagcgtattcaacaaggggctgaaggatgcccagaaggtac 4560
CTCAGTTTTACCGAGAGGAGTTCGCATAAGTTGTTCCCCGACTTCTACGGGTCTTCCATG

4561 cccattgtatgggatctgatctggggcctcggtgcacatgctttacatgtgtttagtcga 4620
GGGTAACATAACCTAGACTAGACCCCGGAGCCACGTGTACGAAATGTACACAAATCAGCT

4621 ggttaaaaaaacgtctaggccccccgaaccacggggacgtggttttcctttgaaaaacac 4680
CCAATTTTTTTGAGATCCGGGGGGCTTGGTGCCCTGCACAAAAGGAACTTTTTTGTG

PuroR marker(4696,5292)>>>

4681 gataataccatggccaccgagtacaagcccacgggtgcgcctcgccaccgcgacgacgtc 4740
CTATTATGGTACCGGTGGCTCATGTTCCGGGTGCCACGCGGAGCGGTGGGCGCTGCTGCAG

4741 cccggggcgtacgcaccctcgccgcccgttcgcccactacccgcccacgcccacacc 4800
GGGGCCCGCATGCGTGGGAGCGGCGGCGCAAGCGGCTGATGGGGCGGTGCGCGGTGTGG

4801 gtcgaccggaccgcccacatcgagcgggtcaccgagctgcaagaactcttctcaccgccc 4860
CAGCTGGGCCTGGCGGTGTAGCTCGCCCAGTGGCTCGACGTTCTTGAGAAGGAGTGCAGC

4861 gtcgggctcgacatcggcaaggtgtgggtcgccgacgacggcgccgcggtggcggtctgg 4920
CAGCCCAGCTGTAGCCGTTCCACACCCAGCGCCTGCTGCCGCGGCGCCACCGCCAGACC

4921 accacgcccggagagcgtcgaagcggggcggtgttcgcccagatcggctcgcgcatggcc 4980
TGGTGGCGCCTCTCGCAGCTTCGCCCCGCCACAAGCGGCTCTAGCCGAGCGCGTACCGG

4981 gagttgagcgggtcccggtgcccgcgcagcaacagatggaaggcctcctggcgcccac 5040
CTCAACTCGCCAAGGGCCGACCGGCGCGTGTGTCTACCTTCCGGAGGACCGCGGCGTG

5041 cggcccaaggagcccgcgtggttccctggccaccgtcggcgtctcgcccaccaccagggc 5100
GCCGGGTTCTCGGGCGCACCAAGGACCGGTGGCAGCCGCAGAGCGGGCTGGTGGTCCCG

5101 aagggctctgggcagcgcctcgtgctccccgagtggaaggcggccgagcgcgctgggggtg 5160
TTCCCAGACCCGTCGCGGCAGCACGAGGGGCTCACCTCCGCCGGCTCGCGCAGCCCCAC

5161 cccgccttctggagacctccgcgccccgcaacctcccccttctacgagcggctcggttc 5220
GGGCGGAAGGACCTCTGGAGGCGCGGGGCGTTGGAGGGGAAGATGCTCGCCGAGCCGAAG

5221 accgtcaccgcccacgtcgaggtgcccgaaggaccgcccacctgggtgcatgaccgcaag 5280
TGGCAGTGGCGGCTGCAGCTCCACGGGCTTCTGGCGCGTGGACCACGTACTGGGCGTTC

```

                    5'mir30(vector_portion) reg(5296,5390)>>>
                    |
5281  cccggtgcctgagtttgtttgaatgaggcttcagtactttacagaatcgttgcctgcaca 5340
      GGGCCACGGACTCAAACAAACTTACTCCGAAGTCATGAAATGTCTTAGCAACGGACGTGT

                                                    XhoI
                                                    |
                    HpaI
                    |
5'mir30(inserted_with_hairpin) reg(5391,5423)>>>
5341  tcttggaaacacttgctgggattacttcttcaggttaaccaacagaaggctcgagAAGG 5400
      AGAACCTTTGTGAACGACCCTAATGAAGAAGTCCAATTGGGTTGTCTTCCGAGCTCTTCC

                                                    mir30-loop
reg(5446,5464)>>>
5401  TATATTGCTGTTGACAGTGAGCGACCTCCACCCTCACTCTGCCATTAGTGAAGCCACAGA 5460
      ATATAACGACAACCTGTCACTCGCTGGAGGTGGGAGTGAGACGGTAATCACTTCGGTGTCT

                                                    3'mir30(vector_portion)
reg(5507,5614)>>>
                    |
                    3'mir30(inserted_with_hairpin)
reg(5487,5506)>>>
5461  TGTAATGGCAGAGTGAGGGTGGAGGGTGCCTACTGCCTCGgaattcaaggggctacttta 5520
      ACATTACCGTCTCACTCCCACCTCCCACGGATGACGGAGCCTTAAGTTCCCCGATGAAAT

5521  ggagcaattatcttgtttactaaaactgaataccttgctatctctttgatacatttttac 5580
      CCTCGTTAATAGAACAAATGATTTTGACTTATGGAACGATAGAGAAACTATGTAAAAATG

                                                    5' Common-Barcode-Flank
other(5626,5646)>>>
5581  aaagctgaattaaatgggtataaattaaatcacttttttcaattggaagactaatgcggc 5640
      TTTCGACTTAATTTTACCATATTTAATTTAGTGAAAAAAGTTAACCTTCTGATTACGCCG

5641  cggccattactccgtctcgtgtcttggtgcatatgtctgctggtttggttgatggtggtt 5700
      GCCGGTAATGAGGCAGAGCACAGAACAACGTATACAGACGACCAAACAACTACAACAAA

                    T7 prom(5710,5729)<<<                    WPRE (HIV)
reg(5749,6337)>>>
                    |                    |
                    3' Common-Barcode-Flank other(5707,5729)>>>
                    |                    |
5701  gcgggcgggccctatagtgagtcgtattacctaggacgcgtctggaacaatcaacctctg 5760
      CGCCCCCGGGATATCACTCAGCATAATGGATCCTGCGCAGACCTTGTTAGTTGGAGAC

5761  gattacaaaatttgtgaaagattgactggattcttaactatgttgctccttttacgcta 5820
      CTAATGTTTTAAACACTTTCTAACTGACCATAAGAATTGATACAACGAGGAAAATGCGAT

5821  tgtggatagctgctttaatgcctttgtatcatgctattgcttcccgatggctttcatt 5880
      ACACCTATGCGACGAAATTACGGAAACATAGTACGATAACGAAGGCATACCGAAAGTAA

5881  ttctcctccttgataaatacctgggtgctgtctctttatgaggagttgtggcccgttgtc 5940

```


AAGAGGAGGAACATATTTAGGACCAACGACAGAGAAATACTCCTCAACACCGGGCAACAG

5941 aggcaacgtggcgtgggtgtgactgtggttctgacgcaacccccactgggtggggcatt 6000
TCCGTTGCACCGCACCACACGTGACACAAACGACTGCGTTGGGGGTGACCAACCCCGTAA

6001 gccaccacctgtcagctcctttccgggactttcgctttccccctccctattgccacggcg 6060
CGGTGGTGGACAGTCGAGGAAAGGCCCTGAAAGCGAAAGGGGGAGGGATAACGGTGCCGC

6061 gaactcatcgccgctgccttgccttgcctgctggacaggggctcggctgttgggactgac 6120
CTTGAGTAGCGGCGGACGGAACGGGCGACGACCTGTCCCGAGCCGACAACCCGTGACTG

6121 aattccggtggtgttgcggggaagctgacgtcctttccatggctgctgcctgtgttggc 6180
TTAAGGCACCACAACAGCCCTTCGACTGCAGGAAAGGTACCGACGAGCGGACACAACGG

6181 acctggattctgcgcgggacgtccttctgctacgtcccttcggccctcaatccagcggac 6240
TGGACCTAAGACGCGCCCTGCAGGAAGACGATGCAGGAAGCCGGGAGTTAGGTCGCCTG

6241 cttccttcccgggctgctgcccggctctgcccctcttccgcgtcttgccttgcct 6300
GAAGGAAGGGCGCCGGACGACGGCCGAGACGCCGAGAAGGCGCAGAAGCGGAAGCGGGA

6301 cagacgagtcggatctcccttggggcgcctccccgcctggaattaattctgcagtcgag 6360
GTCTGCTCAGCCTAGAGGGAAACCCGGCGGAGGGGCGGACCTTAATTAAGACGTGAGCTC

6361 acctagaaaaacatggagcaatcacaagtagcaatacagcagctaccaatgctgattgtg 6420
TGGATCTTTTTGTACCTCGTTAGTGTTCATCGTTATGTCGTCGATGGTTACGACTAACAC

6421 cctggctagaagcacaagaggaggagggtgggttttccagtcacacctcaggtacctt 6480
GGACCGATCTTCGTGTTCTCCTCCTCCACCCAAAAGGTCAGTGTGGAGTCCATGGAA

6481 taagaccaatgacttacaagcagctgtagatcttagccactttttaaagaaaagaggg 6540
ATTCTGTTACTGAATGTTCCGTCGACATCTAGAATCGGTGAAAAATTTTCTTTTCTCCC

3'SIN-LTR(Lenti) other(6544,6779)>>>
|
delta_U3 reg(6544,6596)>>> R(HIV-LTR)
reg(6599,6693)>>>
|

6541 gactggaagggctaattcactcccaacgaagacaagatctgctttttgcttgtactgggt 6600
CTGACCTTCCCGATTAAGTGAGGGTTGCTTCTGTTCTAGACGAAAACGAACATGACCCA

6601 ctctctgggttagaccagatctgagcctgggagctctctgggtaactagggaaaccactgc 6660
GAGAGACCAATCTGGTCTAGACTCGGACCCTCGAGAGACCGATTGATCCCTTGGGTGACG

U5(HIV-LTR) reg(6694,6778)>>>
|

6661 ttaagcctcaataaagcttgccttgagtgttcaagtagtgtgtgcccgtctgttgtgtg 6720
AATTCGGAGTTATTTTGAACGGAACCTCACGAAGTTCATCACACACGGGCAGACAACACAC

6721 actctggtaactagagatccctcagacccttttagtcagtggtgaaaatctctagcagta 6780
TGAGACCATTGATCTCTAGGGAGTCTGGGAAAATCAGTCACACCTTTTAGAGATCGTCAT

6781 gtagttcatgtcatcttattattcagttttataacttgcaaagaaatgaatatcagaga 6840
CATCAAGTACAGTAGAATAATAAGTCATAAATATTGAACGTTTCTTTACTTATAGTCTCT

BGH-polyA reg(6868,7091)>>>
|

```

6841  gtgagaggccttgacattgttttaaacccgctgatcagcctcgactgtgccttctagttgc 6900
      CACTCTCCGGAAGTGTAAACAAATTTGGGCGACTAGTCGGAGCTGACACGGAAGATCAACG

6901  cagccatctgttgtttgcccctccccgctgccttccttgaccctggaaggtgccactccc 6960
      GTCGGTAGACAACAAACGGGGAGGGGGCACGGAAGGAACTGGGACCTTCCACGGTGAGGG

6961  actgtcctttcctaataaaaatgaggaaattgcatcgattgtctgagtaggtgtcattct 7020
      TGACAGGAAAGGATTATTTTACTCCTTTAACGTAGCGTAAACAGACTCATCCACAGTAAGA

7021  attctggggggtgggggtggggcaggacagcaagggggaggattgggaagacaatagcagg 7080
      TAAGACCCCCACCCACCCCGTCTGTCTGTTCCCCCTCCTAACCTTCTGTTATCGTCC

7081  catgctggggatgcggtgggctctatggcttctgaggcggaaagaaccagctggggctct 7140
      GTACGACCCCTACGCCACCCGAGATACCGAAGACTCCGCCTTTCTTGGTTCGACCCCGAGA

      f1 origin(7171,7477)>>>
      |
7141  aggggggatccccacgcgcctgtagcggcgcattaagcgcggcgggtgtggtggttacg 7200
      TCCCCCATAGGGGTGCGCGGGACATCGCCGCGTAATTCGCGCCGCCACACCACCAATGC

7201  cgcagcgtgaccgctacacttgccagcgccttagcgcggcgtcctttcgctttcttcct 7260
      GCGTCGCACTGGCGATGTGAACGGTTCGCGGGATCGCGGGCGAGGAAAGCGAAAGAAGGGA

7261  tcctttctcgccacgttcgcccgtttccccgtcaagctctaaatcgggggctcccttta 7320
      AGGAAAGAGCGGTGCAAGCGGCCGAAAGGGGCAGTTCGAGATTTAGCCCCGAGGGAAT

7321  ggggtccgatttagtgctttacggcacctcgacccccaaaaaacttgattagggtgatggt 7380
      CCCAAGGCTAAATCACGAAATGCCGTGGAGCTGGGGTTTTTTGAACTAATCCCCTACCA

7381  tcacgtagtgggccatcgccctgatagacgggtttttcgcccttgacgttgaggatccag 7440
      AGTGCATCACCCGGTAGCGGGACTATCTGCCAAAAAGCGGGAAACTGCAACCTCAGGTGC

7441  ttctttaatagtgactcttggttccaaactggaacaacactcaaccctatctcggctctat 7500
      AAGAAATTATCACCTGAGAACAAGTTTTGACCTTGTGTGAGTTGGGATAGAGCCAGATA

7501  tcttttgatttataagggattttgcccatttcggcctattgggttaaaaaatgagctgatt 7560
      AGAAAATAAATATTCCCTAAAACGGCTAAAGCCGGATAACCAATTTTTTACTCGACTAA

```

SV40-Enhancer (DTS)

other (7609, 7680) >>>

|
SV40 prom(7609, 7930) >>>

```

7561  taacaaaaatttaacgcgaattaattctgtggaatgtgtgtcagttagggtgtggaaagt 7620
      ATTGTTTTTAAATTGCGCTTAATTAAGACACCTTACACACAGTCAATCCCACACCTTTCA

7621  ccccaggctccccagcaggcagaagtatgcaaagcatgcatctcaattagtcagcaacca 7680
      GGGGTCCGAGGGGTGTCGGTCTTCATACGTTTCGTACGTAGAGTTAATCAGTCGTTGGT

```

SV40-Enhancer (DTS) other (7681, 7752) >>>

```

7681  ggtgtggaaagtccccaggctccccagcaggcagaagtatgcaaagcatgcatctcaatt 7740
      CCACACCTTTAGGGGTCCGAGGGGTGTCGGTCTTCATACGTTTCGTACGTAGAGTTAA

```

SV40 origin(7776, 7853) >>>

7741 agtcagcaaccatagtcgccccctaaactccgccccatccccccccctaaactccgccccagtt 7800
TCAGTCGTTGGTATCAGGGCGGGGATTGAGGCGGGTAGGGCGGGGATTGAGGCGGGTCAA

7801 ccgcccattctccgccccatggctgactaattttttttatgtatgcagaggccgaggccg 7860
GGCGGGTAAGAGGCGGGGTACCGACTGATTAATAAATAAATAACGTCTCCGGCTCCGGC

7861 cctctgcctctgagctattccagaagtagtgaggaggcttttttggaggcctaggctttt 7920
GGAGACGGAGACTCGATAAGGTCTTCATCACTCCTCCGAAAAACCTCCGGATCCGAAAA

HygroR

marker(7979,8996)>>>

7921 gcaaaaagctcccgggagcttgatatccattttcggatctgatcagcacgtgatgaaaa 7980
CGTTTTTCGAGGGCCCTCGAACATATAGGTAAGCCTAGACTAGTCGTGCACTACTTTT

7981 agcctgaactcaccgcgacgtctgtcgagaagtttctgatcgaaaagttcgacagcgtct 8040
TCGGACTTGAGTGGCGCTGCAGACAGCTCTTCAAAGACTAGCTTTTTCAAGCTGTGCAGA

8041 ccgacctgatgcagctctcggagggcgaagaatctcgtgctttcagcttcgatgtaggag 8100
GGCTGGACTACGTGAGAGCCTCCCGCTTCTTAGAGCACGAAAGTCGAAGCTACATCCTC

8101 ggcgtggatatgtcctgccccgtaaatagctgccccgatggtttctacaaagatcgttatg 8160
CCGCACCTATACAGGACGCCATTTATCGACGCGGCTACCAAAGATGTTTTCTAGCAATAC

8161 tttatcggcactttgcatcgcccgcgctcccgattccggaagtgcttgacattggggaaat 8220
AAATAGCCGTGAAACGTAGCCGGCGCGAGGGCTAAGGCCTTACGAAGTGTAAACCCCTTA

8221 tcagcgagagcctgacctattgcatctcccgcgctgcacaggggtgtcacgttgcaagacc 8280
AGTCGCTCTCGGACTGGATAACGTAGAGGGCGGCACGTGTCCACAGTGCAACGTTCTGG

8281 tgcctgaaaccgaactgccccgctgttctgcagccggtcgcggaggccatggatgcgatcg 8340
ACGGACTTTGGCTTGACGGGCGACAAGACGTGCGCCAGCGCCTCCGGTACCTACGCTAGC

8341 ctgccccgatcttagccagacgagcgggttcggccattcggaccgcaaggaatcggtc 8400
GACGCCGGCTAGAATCGGTCTGCTCGCCCAAGCCGGTAAGCCTGGCGTTCCTTAGCCAG

8401 aatacactacatggcgtgatttcatatgcgcgattgctgatccccatgtgtatcactggc 8460
TTATGTGATGTACCGCACTAAAGTATACGCGCTAACGACTAGGGGTACACATAGTGACCG

8461 aaactgtgatggacgacaccgtcagtgcgctccgtcgcgcaggctctcgatgagctgatgc 8520
TTTGACACTACCTGCTGTGGCAGTCACGCAGGCAGCGCTCCGAGAGCTACTCGACTACG

8521 tttgggcccaggactgccccgaagtccggcacctcgtgcacgcggatttcggctccaaca 8580
AAACCCGGCTCCTGACGGGGCTTACGGCCGTGGAGCACGTGCGCCTAAAGCCGAGGTTGT

8581 atgtcctgacggacaatggccgcataacagcggctcattgactggagcggagcgtatgtcg 8640
TACAGGACTGCCTGTTACCGGCGTATTGTGCGCCAGTAACTGACCTCGCTCCGCTACAAGC

8641 gggattcccaatacagaggtcgccaacatcttcttctggaggccgtgggttgcttgatgg 8700
CCCTAAGGGTTATGCTCCAGCGGTTGTAGAAGAAGACCTCCGGCACCAACCGAACATACC

8701 agcagcagacgcgctacttcgagcggaggcatccggagcttgaggatcgcccgggctcc 8760
TCGTCGCTGCGCGATGAAGCTCGCCTCCGTAGGCCTCGAACGTCCTAGCGGCGCCGAGG

8761 gggcgatatatgctccgcattggctcttgaccaactctatcagagcttggttgacggcaatt 8820
CCCGCATATACGAGGCGTAACCGAAGTGGTTGAGATAGTCTCGAACCAACTGCCGTTAA

8821 tcgatgatgcagcttggggcgcagggtcgatgacgacgcaatcgccgatccggagccggga 8880
AGCTACTACGTGCAACCCGCGTCCCAGCTACGCTGCGTTAGCAGGCTAGGCCTCGGCCCT

8881 ctgtcggggtacacaaaatcgcccgcagaagcgcggccgtctggaccgatggctgtgtag 8940
GACAGCCCGCATGTGTTTTAGCGGGCGTCTTCGCGCCGGCAGACCTGGCTACCGACACATC

8941 aagtactcgccgatagtggaaaccgacgccccagcactcgtccgagggcaaaggaatagc 9000
TTCATGAGCGGCTATCACCTTTGGCTGCGGGGTCGTGAGCAGGCTCCCGTTTCCTTATCG

9001 acgtgctacgagatttcgattccaccgccccttctatgaaagggtgggcttcggaatcg 9060
TGCACGATGCTCTAAAGCTAAGGTGGCGGCGGAAGATACTTTCCAACCCGAAGCCTTAGC

9061 ttttcggggacgcccggctggatgatcctccagcgcggggatctcatgctggagttcttcg 9120
AAAAGGCCCTGCGGCCGACCTACTAGGAGGTGCGCCCCCTAGAGTACGACCTCAAGAAGC

SV40-polyA-signal

reg(9160,9194)>>>

9121 cccaccccaacttgtttattgcagcttataatggttacaaataaagcaatagcatcacia 9180
GGGTGGGGTTGAACAAATAACGTGCAATATTACCAATGTTTATTTTCGTTATCGTAGTGTT

9181 atttcaciaataaagcatttttttactgcattctagttgtggtttgtccaaactcatca 9240
TAAAGTGTTTATTTTCGTAAAAAAGTGACGTAAGATCAACACCAAACAGGTTTGAGTAGT

9241 atgtatcttatcatgtctgtataaccgtcgacctctagctagagcttggcgtaatcatggt 9300
TACATAGAATAGTACAGACATATGGCAGCTGGAGATCGATCTCGAACCGCATTAGTACCA

lac prom(9342,9425)<<<

9301 catagctgtttcctgtgtgaaattggtatccgctcacaattccacacaacatacagaccg 9360
GTATCGACAAAGGACACACTTTAACAATAGGCGAGTGTTAAGGTGTGTTGTATGCTCGGC

9361 gaagcataaagtgtaaagcctgggggtgcctaataagtgagtgagtaactcacattaattgct 9420
CTTCGTATTTTACATTTTCGGACCCACGATTACTCACTCGATTGAGTGTAATTAACGCA

9421 tgcgctcactgcccgctttccagtcgggaaacctgtcgtgccagctgcattaatgaatcg 9480
ACGCGAGTGACGGGCGAAAGGTCAGCCCTTTGGACAGCACGGTCGACGTAATTACTTAGC

9481 gccaacgcgcggggagaggcgggtttgctgattggggcgtcttccgcttcctcgtcactg 9540
CGGTTGCGCGCCCTCTCCGCCAAACGCATAACCCGCGAGAAGGCGAAGGAGCGAGTGAC

9541 actcgtcgcgctcggctcgttcggctgcggcgagcggatcagctcactcaaaggcggtaa 9600
TGAGCGACGCGAGCCAGCAAGCCGACGCCGCTCGCCATAGTCGAGTGAGTTTCCGCCATT

9601 tacggttatccacagaatcaggggataacgcaggaagaacatgtgagcaaaaggccagc 9660
ATGCCAATAGGTGTCTTAGTCCCCTATTGCGTCCTTTCTTGTACACTCGTTTTCCGGTCG

pUC origin(9686,10305)<<<

9661 aaaaggccaggaaccgtaaaaaggcgcggttgcgctttttccataggctccgcccc 9720
TTTTCCGGTCCTTGGCATTTTTCCGGCGCAACGACCGCAAAAAGGTATCCGAGGCGGGG

9721 ctgacgagcatcacaaaaatcgacgctcaagtgcagaggtggcgaaacccgacaggactat 9780
GACTGCTCGTAGTGTTTTTAGCTGCGAGTTCAGTCTCCACCGCTTTGGGCTGTCCTGATA

9781 aaagataaccaggcggtttccccctggaagctccctcgtgcgctctcctgttccgaccctgc 9840
TTTCTATGGTCCGCAAAGGGGGACCTTCGAGGGAGCACGCGAGAGGACAAGGCTGGGACG

9841 cgcttaccggatacctgtccgcctttctcccttcgggaagcgtggcgctttctcatagct 9900
GCGAATGGCCTATGGACAGGCGGAAAGAGGGAAGCCCTTCGCACCGCGAAAGAGTATCGA

9901 cacgctgtaggtatctcagttcgggtgtaggtcgctcctccaagctgggctgtgtgcacg 9960
GTGCGACATCCATAGAGTCAAGCCACATCCAGCAAGCGAGGTTTCGACCCGACACACGTGC

9961 aaccccccggttcagccccgaccgctgcgcccttatccggtaactatcgtcttgagtccaacc 10020
TTGGGGGGCAAGTCGGGCTGGCGACGCGGAATAGGCCATTGATAGCAGAACTCAGGTTGG

10021 cggtaagacacgacttatcgccactggcagcagccactggtaacaggattagcagagcga 10080
GCCATTCTGTGCTGAATAGCGGTGACCGTCGTCGGTGACCATTGTCCTAATCGTCTCGCT

10081 ggtatgtaggcgggtgctacagagttcttgaagtgggtggcctaactacggctacactagaa 10140
CCATACATCCGCCACGATGTCTCAAGAACTTCACCACCGGATTGATGCCGATGTGATCTT

10141 gaacagtatattgggtatctgcgctctgctgaagccagttaccttcggaaaaagagttggta 10200
CTTGTCTATAAACCATAGACGCGAGACGACTTCGGTCAATGGAAGCCTTTTTCTCAACCAT

10201 gctcttgatccggcaaacacaccaccgctggtagcgggtgggtttttttgtttgcaagcagc 10260
CGAGAACTAGGCCGTTTTGTTTGGTGGCGACCATCGCCACCAAAAAACAAACGTTTCGTCG

10261 agattacgcgcagaaaaaaaggatctcaagaagatcctttgatcttttctacggggctcg 10320
TCTAATGCGCGTCTTTTTTCTTAGAGTTCTTCTAGGAACTAGAAAAGATGCCCCAGAC

10321 acgctcagtggaacgaaaaactcacgttaagggattttgggtcatgagattatcaaaaagga 10380
TGCGAGTCACTTTGCTTTTGGAGTGCAATTCCTAAAACAGTACTCTAATAGTTTTTCT

10381 tcttcacctagatccttttaaatataaaatgaagttttaaatcaatctaaagtatatatg 10440
AGAAGTGGATCTAGGAAAATTTAATTTTTACTTCAAATTTAGTTAGATTTTCATATATAC

AmpR marker(10460,11320)<<<
|

10441 agtaaaacttgggtctgacagttaccaatgcttaatcagtgaggcacctatctcagcgatct 10500
TCATTTGAACCAGACTGTCAATGGTTACGAATTAGTCACTCCGTGGATAGAGTCGCTAGA

10501 gtctatctcggttcatccatagttgcctgactccccgctcgtgtagataactacgatacggg 10560
CAGATAAAGCAAGTAGGTATCAACGGACTGAGGGGCAGCACATCTATTGATGCTATGCC

10561 agggcttaccatctggccccagtgctgcaatgataccgcgagaccacgctcaccggctc 10620
TCCCGAATGGTAGACCGGGGTCACGACGTTACTATGGCGCTCTGGGTGCGAGTGGCCGAG

10621 cagatttatcagcaataaaccagccagccggaagggccgagcgcagaagtggctcctgcaa 10680
GTCTAAATAGTCGTTATTTGGTTCGGTTCGGCCTTCCCGGCTCGCGTCTTCACCAGGACGTT

10681 ctttatccgcctccatccagctctattaattggtgcccgggaagctagagtaagttagttcgc 10740
GAAATAGGCGGAGGTAGGTCAGATAATTAACAACGGCCCTTCGATCTCATTCAATCAAGCG

10741 cagttaatagtttgcgcaacggttggtgcccattgctacagggcatcgtgggtgtcacgctcgt 10800
GTCAATTATCAAACGCGTTGCAACAACGGTAACGATGTCCGTAGCACCACAGTGCAGCA

10801 cgtttggtatggcttcattcagctccgggtcccaacgatcaaggcgagttacatgatccc 10860
GCAAACCATAACGAAGTAAGTCGAGGCCAAGGGTTGCTAGTTCGGCTCAATGTACTAGGG

10861 ccatgttgtgcaaaaaagcggttagctccttcgggcctccgatcgttgtcagaagtaagt 10920
GGTACAACACGTTTTTTTCGCCAATCGAGGAAGCCAGGAGGCTAGCAACAGTCTTCATTCA

10921 tggccgcagtggtatcactcatgggtatggcagcactgcataattctcttactgtcatgc 10980
ACCGGCGTCACAATAGTGAGTACCAATACCGTCGTGACGTATTAAGAGAATGACAGTACG

10981 catccgtaagatgcttttctgtgactggtgagtactcaaccaagtcattctgagaatagt 11040
GTAGGCATTCTACGAAAAGACACTGACCACTCATGAGTTGGTTCAGTAAGACTCTTATCA

11041 gtatgccccgaccgagttgctcttgccccggcgtcaatacgggataataccgcgccacata 11100
CATACGCCGCTGGCTCAACGAGAACGGGCCGCAGTTATGCCCTATTATGGCGCGGTGTAT

11101 gcagaactttaaaagtgctcatcattggaaaacgttcttcggggcgaaaactctcaagga 11160
CGTCTTGAAATTTTTCACGAGTAGTAACCTTTTGCAGAAGCCCCGCTTTTGAGAGTTCCT

11161 tcttaccgctggtgagatccagttcgatgtaaccactcgtgcacccaactgatcttcag 11220
AGAATGGCGACAACCTTAGGTCAAGCTACATTGGGTGAGCACGTGGGTTGACTAGAAGTC

11221 catcttttactttcaccagcgtttctgggtgagcaaaaacaggaaggcaaaatgccgcaa 11280
GTAGAAAATGAAAGTGGTTCGCAAAGACCCACTCGTTTTTGTCTTCCGTTTTACGGCGTT

11281 aaaaggggaataagggcgacacggaaatggtgaatactcatactcttcctttttcaatatt 11340
TTTTCCCTTATTCGCTGTGCCTTTACAACCTTATGAGTATGAGAAGGAAAAAGTTATAA

Bacterial-Promoter prom(11362,11400)<<<
|

11341 attgaagcatttatcagggttattgtctcatgagcggatacatatttgaatgtatttaga 11400
TAACTTCGTAAATAGTCCCAATAACAGAGTACTCGCCTATGTATAAACTTACATAAATCT

11401 aaaataaacaaataggggttccgcgcacatttccccgaaaagtgccacctgacgtcgagc 11460
TTTTATTTGTTTATCCCCAAGGCGCGTGTAAAGGGGCTTTTTCACGGTGGACTGCAGCTGC

SV40-polyA-signal

reg(11504,11538)>>>

|

11461 gatcgggagatcaacttgtttattgcagcttataatgggttaciaaataaagcaatagcatc 11520
CTAGCCCTCTAGTTGAACAAATAACGTCGAATATTACCAATGTTTATTTTCGTTATCGTAG

11521 acaaatttcacaaataaagcatttttttctactgcattctagttgtggtttgcctcaaactc 11580
TGTTTAAAGTGTATTTTCGTAAAAAAGTGACGTAAGATCAACACCAAACAGGTTTGAG

11581 atcaatgtatcttatcatgtctggatcaactggataactcaagctaaccaaataatcatccc 11640
TAGTTACATAGAATAGTACAGACCTAGTTGACCTATTGAGTTCGATTGGTTTTAGTAGGG

11641 aaacttcccaccccataccctattaccactgccaattacctgtggtttctttactctaa 11700
TTTGAAGGGTGGGGTATGGGATAATGGTGACGGTTAATGGACACCAAAGTAAATGAGATT

11701 acctgtgattcctctgaattattttcatttttaagaaattgtatttgtaaatatgtact 11760
TGGACACTAAGGAGACTTAATAAAAAGTAAAATTTCTTTAACATAAAACAATTTATACATGA

11761 acaaacttagtagt 11774
TGTTTGAATCATCA

Restriction analysis of pGIPZ lentiviral vector

AhdI (GACnn_n'nnGTC) [Eam1105I,AspEI,DriI,EclHKI]
Cuts 1 time.
Cuts at position 10533.
Fragment sizes 10533, 1241.

AleI (CACnn'nnGTG) [OliI]
Cuts 1 time.
Cuts at position 1577.
Fragment sizes 1577, 10197.

AloI (GAACnnnnnnTCCnnnnnnn_nnnnn')
Cuts 1 time.
Cuts at position 7423.
Fragment sizes 7423, 4351.

AloI (GGAnnnnnnGTTcnnnnnnn_nnnnn')
Cuts 1 time.
Cuts at position 7455.
Fragment sizes 7455, 4319.

AsiSI (GCG_AT'CGC) [SgfI]
Cuts 1 time.
Cuts at position 8338.
Fragment sizes 8338, 3436.

BbvCI (CC'TCA_GC)
Cuts 1 time.
Cuts at position 1424.
Fragment sizes 1424, 10350.

BlpI (GC'TnA_GC) [Bpu1102I,Bsp1720I,CelII]
Cuts 1 time.
Cuts at position 3564.
Fragment sizes 3564, 8210.

Bpu10I (CC'TnA_GC)
Cuts 1 time.
Cuts at position 1424.
Fragment sizes 1424, 10350.

BsaBI (GATnn'nnATC) [Bse8I,BseJI,MamI]
[dam methylated]
Cuts 1 time.
Cuts at position [3853].
Fragment sizes 3853, 7921.

BsiWI (C'GTAC_G) [Pfl23II,PspLI,SunI]

Cuts 1 time.
Cuts at position 4749.
Fragment sizes 4749, 7025.

BsrGI (T'GTAC_A) [Bsp1407I,BstAUI,SspBI]
Cuts 1 time.
Cuts at position 4089.
Fragment sizes 4089, 7685.

BstEII (G'GTnAC_C) [BstPI,Eco91I,EcoO65I,PspEI]
Cuts 1 time.
Cuts at position 4827.
Fragment sizes 4827, 6947.

BstZ17I (GTA'TAC) [BssNAI,Bst1107I]
Cuts 1 time.
Cuts at position 9261.
Fragment sizes 9261, 2513.

Bsu36I (CC'TnA_GG) [AxyI,Bse21I,Eco81I]
Cuts 1 time.
Cuts at position 6469.
Fragment sizes 6469, 5305.

CspCI (CAAnnnnnGTGnnnnnnnnnn_nn')
Cuts 1 time.
Cuts at position 3141.
Fragment sizes 3141, 8633.

CspCI (CCACnnnnnTTGnnnnnnnnnn_nn')
Cuts 1 time.
Cuts at position 3106.
Fragment sizes 3106, 8668.

EcoNI (CCTnn'n_nnAGG) [BstENI,XagI]
Cuts 1 time.
Cuts at position 1170.
Fragment sizes 1170, 10604.

FspI (TGC'GCA) [Acc16I,AviII,NsbI]
Cuts 1 time.
Cuts at position 10755.
Fragment sizes 10755, 1019.

HpaI (GTT'AAC) [KspAI]
Cuts 1 time.
Cuts at position 5376.
Fragment sizes 5376, 6398.

MluI (A'CGCG_T)
Cuts 1 time.

Cuts at position 5736.
Fragment sizes 5736, 6038.

NotI (GC'GGCC_GC) [CciNI]

Cuts 1 time.
Cuts at position 4100.
Fragment sizes 4100, 7674.

NruI (TCG'CGA) [Bsp68I]

[dam methylated]
Cuts 1 time.
Cuts at position [833].
Fragment sizes 833, 10941.

PmeI (GTTT'AAAC) [MssI]

Cuts 1 time.
Cuts at position 6862.
Fragment sizes 6862, 4912.

PpuMI (rG'GwC_Cy) [PpuXI,Psp5II,PspPPI]

[dcm methylated]
Cuts 1 time.
Cuts at position 1934.
Fragment sizes 1934, 9840.

PshAI (GACnn'nnGTC) [BoxI,BstPAI]

Cuts 1 time.
Cuts at position 8001.
Fragment sizes 8001, 3773.

SanDI (GG'GwC_CC)

Cuts 1 time.
Cuts at position 1934.
Fragment sizes 1934, 9840.

SfiI (GGCCn_nnn'nGGCC)

[dcm methylated]
Cuts 1 time.
Cuts at position 2621.
Fragment sizes 2621, 9153.

SgrAI (Cr'CCGG_yG)

Cuts 1 time.
Cuts at position 2500.
Fragment sizes 2500, 9274.

SnaBI (TAC'GTA) [BstSNI,Eco105I]

Cuts 1 time.
Cuts at position 3070.
Fragment sizes 3070, 8704.

SspI (AAT'ATT)
Cuts 1 time.
Cuts at position 11337.
Fragment sizes 11337, 437.

XbaI (T'CTAG_A)
[dam methylated]
Cuts 1 time.
Cuts at position 2707.
Fragment sizes 2707, 9067.

XhoI (C'TCGA_G) [BssHI, PaeR7I, Sfr274I, SlaI, StrI, TliI]
Cuts 1 time.
Cuts at position 5391.
Fragment sizes 5391, 6383.

References

- Bartel, DP (2004)** "MicroRNAs: genomics, biogenesis, mechanism and function" *Cell*, 116, 281–297.
- Boden, D et al (2004)** "Enhanced Gene silencing of HIV-1 specific siRNA using microRNA designed hairpins" *Nucleic Acids Research* 32:3, 1154–58.
- Cullen, BR et al (2005)** "RNAi the natural way" *Nature Genetics* 37:11, 1163-1165.
- Cullen, BR (2004)** "Transcription and processing of human microRNA precursors" *Mol. Cell* 16, 861–865.
- Cleary, M et al (2004)** "Production of complex nucleic acid libraries using highly parallel in situ oligonucleotide synthesis" *Nature Methods* 1:3, 241-248.
- Chendrimada, TP et al (2005)** "TRBP recruits the Dicer complex to Ago2 for microRNA processing and gene silencing" *Nature* 436, 740-744.
- Elbashir, SM et al (2001)** "Duplexes of 21-nucleotide RNAs mediate RNA interference in cultured mammalian cells" *Nature*, 411: 494–498.
- Fire, A et al (1998)** "Potent and specific genetic interference by double stranded RNA in *Caenorhabditis elegans*" *Nature* 391:806-811.
- Gregory, RI et al (2004)** "The Microprocessor complex mediates the genesis of microRNAs" *Nature*, 432(7014):235-40.
- Gregory, RI et al (2005)** "Human RISC couples microRNA biogenesis and posttranscriptional gene silencing" *Cell* 18;123 (4):631-40
- Hannon, GJ et al (2003)** "Stable Suppression of Gene Expression by RNAi in Mammalian Cells" *PNAS* 99 (3):1443–1448.
- Hannon, GJ et al (2004)** "Unlocking the potential of the human genome with RNA interference" *Nature* 431:371–378.
- He, L and Hannon, GJ (2004)** "MicroRNAs: Small RNAs with a big role in Gene Regulation" *Nature Genetics Reviews* 5, 522–531.
- Lee, Y et al (2002)** "MicroRNA maturation: stepwise processing and subcellular localization" *EMBO* 21, 17: 4663-4670.
- Lee, Y et al (2003)** "The nuclear RNase III Drosha initiates microRNA processing" *Nature* 425:415–418.
- Nakahara, et al (2004)** "Expanding roles for miRNAs and siRNAs in cell regulation" *Current Opinion in Cell Biology* 16:127–133.

Paddison, PJ et al. (2004) "A resource for large-scale RNAi based screens in mammals"
Nature 428: 427–431.

Silva, JM et al (2005) "Second-generation shRNA libraries covering the mouse and human genomes" *Nat. Genet.* 37:11, 1281–1288.

Stegmeier, F et al (2004) "A lentiviral microRNA-based system for single-copy polymerase II-regulated RNA interference in mammalian cells" *Proc. Natl. Acad. Sci. USA* 102, 13212–13217.

Zeng, Y et al (2003) "MicroRNAs and small interfering RNAs can inhibit mRNA expression by similar mechanisms". *PNAS* 100: 9779–9784.

Zeng, Y et al (2005) "Recognition and cleavage of primary microRNA precursors by the nuclear processing enzyme Drosha" *EMBO J.* 24: 138-148.

Limited Use License:

This product is covered by several patent applications owned by Cold Spring Harbor Laboratory. The purchase of this product conveys to the buyer the limited, non-exclusive, non-transferable right (without the right to resell, repackage, or further sublicense) under these patent rights to perform the RNAi knockdown methods using the RNAi-inducing vectors claimed in those patent applications for research purposes solely in conjunction with this product. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of this product does not include nor carry any right or license to use, develop, or otherwise exploit this product commercially, and no rights are conveyed to the buyer to use the product or components of the product for any other purposes, including without limitation, provision of services to a third party, generation of commercial databases, or clinical diagnostics or therapeutics.

In addition, any commercial organization that purchases or desires to purchase RNAi clones are outside the above research license and should contact Open Biosystems for a research use license.

This product is sold pursuant to a license from CSHL, and CSHL reserves all other rights under these patent rights. For information on purchasing a license to the patent rights for uses other than in conjunction with this product or to use this product for purposes other than research, please contact the CSHL Office of Technology Transfer at, 516-367-8312.