

Pipetting Aid

Lithium battery with no memory effect and up to 12 hours of continuous use

Can be charged while using

Fully autoclavable nose cone and UV resistant housing

Light and fatigue free



Features :

- Lightweight, only 167 grams (including battery)
- A full range of 1-100ml
- Fills 50ml pipet in less than 10 seconds
- Ergonomic design allows longer and fatigue-free pipetting.
- It has a dispense selection of gravity and blow-out with a fully adjustable speed-setting.
- Self-inverted design. It can easily be placed on the bench top with the pipette mounted - without the pipette touching the bench and without risk of backflow or contamination.
- It comes with an exchangeable 0.2 μm PTFE disc filter, acting as shutdown system if overflow occurs, and it has a protective design for venting out corrosive liquids or fumes.
- Express dispensing : depress the dispensing button and pass the built-in point of resistance.
- Silicon made exchangeable pipette holders

Technical Specifications:

Aspirating/Dispensing Speed	> 6ml/sec
Dispense Mode	Blow-out / Gravity Selection
Speed Control Wheel	Continuous
Power Adaptor	100 / 120 / 220 / 240V, AC-DC 3.6V
Pipet Types	Glass or Plastic Pipets(1-100ml)
Filter	Hydrophobic 25mm diameter, 0.2 μm
Autoclavable	Nose cone, pipette holder and filter only
Operating Temperature	+10°C ~ 35°C
Operating Humidity	RH: 0 ~ 85%
Battery	Lithium 600mAh / 3.6V or above
Weight	167g, including a battery
Certification	Complies with CE, Class A (EN61010-1, EN60101-1-2, EN50082-1, EN 55011) RoHS