

An ISO 9001 : 2000 Company



LABTOP[®]
Quality Lab Equipment



The LABTOP[®] Blood Bag Tube Sealer is a compact equipment to seal the Blood Bag pilot tube by radio frequency sealing system. Sealing can be achieved in less than 1.5 seconds per tube by simple operation. Sealed tube can be easily separated by pulling the sealing section by both sides.

BLOOD BAG TUBE SEALER

BLOOD BAG TUBE SEALER



Technical Specification:-

Catalogue No	LTS-3
Oscillation Frequency	40.68 MHz
Radio Frequency Output	15 W
Sealing time (Approx.)	1.5 Sec
Input Voltage	230 V AC
Frequency	50 Hz
Power Consumption	100 V A
Weight Approx	5.5 Kgs

NB: Specifications are subject to change without prior notice

SIE SKYLAB INSTRUMENTS & ENGINEERING PVT. LTD.

M-7 & 8, Anand Mangal Ind. Estate, Survey No.: 70-B,
Opp. IPOL, Waliv Phata, Sativali Road,
Vasai (E), Dist - Thane-401 208, Maharashtra, INDIA.
Tel.: 91-250-2481431, 2481428 Fax : 91-250-2481505
e-mail : byjuskylab@vsnl.net www.skylabinstrument.com

