

Caspase-12 Inhibitor, Z-ATAD-FMK

CATALOG #:	1079-100
LOT #:	_____
STORAGE CONDITIONS:	Store at -20°C.
SHELF LIFE:	6-12 months under proper storage conditions
MOLECULAR WEIGHT:	540
SEQUENCE:	Z-Ala-Thr-Ala-Asp-FMK (FMK, fluoromethyl ketone)
PURITY:	>95% by HPLC analysis.
PACKAGE CONTENTS:	100 µl (2 mM) Z-ATAD-FMK in DMSO

DESCRIPTION:

A synthetic peptide that irreversibly inhibits caspase-12 and related caspase activity. The inhibitor is designed as a methyl ester to facilitate cell permeability.

RECOMMENDED USAGE:

We recommend using 1000 dilutions for inhibiting caspase-12 activity in tissue and cell lysate (Add 1 µl to 1 ml of tissue and cell lysate). The optimal doses may vary for different tissue and cell preparations.

NOTE: For research use only! Not to be used in human.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptotic Cell Isolation Kit

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP & PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Quality Antibodies for Apoptosis and Signal Transduction Molecules