

Caspase-11 Polyclonal Antibody

CATALOG #:	3021-100
AMOUNT:	100 µg
LOT #:	_____
HOST:	Rabbit
IMMUNOGEN:	Synthetic peptide from rat caspase-11 (ID# BV-17)

FORMULATION:

100 µg (0.2 mg/ml) protein A purified rabbit anti-caspase-11 polyclonal antibody in phosphate buffered saline (PBS), pH 7.2, containing 50% glycerol, 1% BSA, 0.02% sodium azide.

SPECIES REACTIVITY: mouse, rat

STORAGE CONDITIONS:

Store at -20°C. For long term storage, aliquot and refreeze at -70°C. Avoid repeated freeze/defrost cycles.

BACKGROUND:

Caspase family of cysteine proteases has been shown to play a key role in apoptosis. Caspase-11 is also synthesized as inactive pro-enzyme that are processed in cells undergoing apoptosis or inflammation. The processed form of caspase-11 consists of large and small subunits which associate to form an active enzyme.

SPECIFICITY:

The antibody detects procaspase-11 (45 kDa), a partially cleaved fragment (30 kDa) and fully cleaved fragment (20 kDa). A 58 kDa unknown identity band was also detected in 3T3 cells.

APPLICATION AND USAGE:

The antibody can be used for Western blotting (1 µg/ml). However, the optimal conditions should be determined individually. Other applications have not been tested.

NOTE: For Research Use Only! Not to be used in human subjects.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies