

Granzyme B Monoclonal Antibody

CATALOG #:	3073-100
LOT #:	_____
HOST:	Mouse
ISOTYPE:	IgG ₁
CLONE:	B18.1
IMMUNOGEN:	Human recombinant granzyme B (a.a. 7-227)
SPECIES REACTIVITY:	Human

FORMULATION:

100 µg (0.5 mg/ml) in PBS containing 1 mg/ml BSA, 1.5 mM sodium azide and 50% glycerol.

STORAGE:

Store at -20°C. For long term storage, aliquot and refreeze at -70°C. Avoid repeated freeze/thaw cycles.

DESCRIPTION:

Granzyme B is a 27 kDa serine protease stored in granules of activated cytotoxic T cells and NK cells. Upon target cell contact, granzyme B is directionally exocytosed and assisted by perforin, enters the target cell. With its unique substrate specificity (granzyme B cleaves after Asp), granzyme B processes and activates various pro-caspases thereby inducing apoptosis in the target cell.

APPLICATIONS:

Western blot, immunoprecipitation, flow cytometry and Immunohistochemistry. However, the optimal conditions should be determined individually.

NOTE: For research use only! Not to be used in human.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Cathepsin and Calpain Assay Kits & Reagents
- Apoptosis Inducers and Set

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione Fluorometric & Colorimetric Assay Kits
- Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP & PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit

Growth Factors and Cytokines (many)