



DATA SHEET

MODEL: LAB/ULT - Ultra Low Temperature Freezers -85 C

Model number:		11 ULT/LAB 11	31 ULT/LAB 31	51 ULT/LAB 51
Ambient temperature	°C	+10° to +30°	+10° to +30°	+10° to +30°
Net volume	litres	130	301	418
Net volume	cu.ft.	4,6	10,6	14,8
Height with 50 mm castors	mm	920	920	920
Height with open lid	mm	1620	1620	1620
Width	mm	720	1300	1700
Insulation thickness	mm	100	100	100
Depth incl. handle and hinges	mm	730	730	730
Inner dimensions (HxDxW)	mm	650x450x 520/320	650x450x 1100/900	650x450x 1500/1300
Power consumption at +25°C	Watts	480	710	810
Fuse	A	10	10	10
Temperature range	°C	-60° to -85°	-60° to -85°	-60° to -85°
Noise level	dBa	<51	<51	<51
Weight	kg	53	78	95
Pull down from +25°C	hours	4	5	5
Handle with lock		+	+	+
Defrost drain		+	+	+
Electric thermometer		+	+	+
Acoustic alarm		+	+	+
Connection to extra alarm available		+	+	+
Castors		+	+	+

Subject to alternations. The technical specifications should be taken as a guidance.