



MONOBODIES

ANTI-HIV-1 p24

Clone 39/5.1.23

ZMC Catalog #: 0801079 **Unit Size:** 100 µg

ZMC Catalog #: 0801126 **Unit Size:** 1 mg

PRODUCT

CHARACTERISTICS

This murine monoclonal antibody reacts with Human Immunodeficiency Virus Type 1 (HIV-1) p24 protein. The antibody

is of the IgG1 subclass and has been purified from serum-free culture supernatant by ammonium sulfate precipitation.

The antibody is supplied as purified immunoglobulin at a concentration of approximately 1 mg/ml in phosphate-buffered

saline. MONOBODIES Anti-HIV-1 p24 (clone 39/5.1.23) exhibits reactivity with HIV-1 infected cultures using indirect

immunofluorescence. It also exhibits reactivity with viral lysates examined using Western blotting. In addition, the

antibody immunoprecipitates p24 protein from radiolabeled infected cell extracts.

RECOMMENDED

USAGE

Antibody dilutions should be prepared using buffers containing suitable protein in order to stabilize antibody activity.

The antibody is used to detect HIV-1 core protein in virus infected cells as well as in viral or infected cell lysates.

Studies on core antigen synthesis and metabolism can be performed using Western blotting or

radioimmunoprecipitation analysis.

STORAGE Aliquot and store at -20°C.

1004160 Rev.: 04/00

All products are for Research Use Only and not for *in vitro* Diagnostic Use

This Package Insert is provided for product evaluation purposes only and is not intended to be used in place of the Package Insert

shipped with the product.

872 Main Street

Buffalo, NY 14202

(800)274-5487 fax (716)882-0959

Trial Version