



MONOBODIES

ANTI-HTLV-I gp 46

Clone 68/4.11.21

ZMC Catalog #: 0801085 **Unit Size:** 100 µg

ZMC Catalog #: 0801118 **Unit Size:** 1 mg

PRODUCT

CHARACTERISTICS

This murine monoclonal antibody reacts with Human T-Lymphotropic Virus Type I (HTLV-I) surface protein, gp46. The epitope is located between amino acids 190-209 of the viral *env* gene. No reactivity with the surface protein of HTLV-II (MoT) has been observed. The antibody is of the IgG1 subclass and has been purified from serum-free culture supernatant by concentration and ammonium sulfate precipitation. The antibody is greater than 85% pure as determined using denaturing polyacrylamide gel electrophoresis. The antibody is supplied as purified immunoglobulin in phosphate-buffered saline, pH 7.2. MONOBODIES Anti-HTLV-I gp46 (clone 68/4.11.21) exhibits reactivity with viral lysates using solid-phase enzyme immunoassay and immunoblotting. It also reacts with HTLV-I infected cell lines using indirect immunofluorescence. The antibody immunoprecipitates gp46 surface protein from radiolabeled cell lysates, as well as viral precursor protein.

RECOMMENDED

USAGE

Antibody dilutions should be prepared using buffers containing suitable protein in order to stabilize antibody activity. The antibody is used to detect HTLV-I surface protein and its precursor in lysates prepared from virions or infected cells. Using immunoblotting and/or radioimmunoprecipitation analysis, the biosynthesis and metabolism of viral surface protein may be studied.

STORAGE Aliquot and store at -20°C.

1004165 Rev.: 04/00

All products are for Research Use Only and not for *in vitro* Diagnostic Use

This Package Insert is provided for product evaluation purposes only and is not intended to be used in place of the Package Insert shipped with the product.

872 Main Street

Buffalo, NY 14202

(800)274-5487 fax (716)882-0959

Trial Version