

PRODUCT: Cell Lysis Buffer

CATALOG #:	1067-100	100 ml
	1067-400	400 ml
LOT #:	_____	
STORAGE:	Store at 4°C	
SHELF LIFE:	One year under proper storage conditions	

DESCRIPTION:

The Cell Lysis Buffer is optimized for lysis of mammalian cells and tissues for caspase activity assays. The buffer is intended to use as a supplemental reagent for researchers who need to lyse large amount of cells or tissues.

For a complete protocol of using the buffer in caspase assays, please refer to BioVision's Caspase Assay Kits (#K105, K106, K112, K113, K114, K115, K116, K117, K118, K119, etc.).

RELATED PRODUCTS:

1. KinaseSTAR™ JNK & Akt Activity Assay Kit
2. Histone Deacetylase Fluorometric & Colorimetric Assay Kits
3. Calpain Activity Assay Kit
4. Apoptotic Cell Isolation Kit
5. MitoCapture™ Mitochondrial Apoptosis Detection Kit
6. Mitochondria/Cytosol Fractionation Kit
7. Nuclear/Cytosol Fractionation Kit
8. Cytochrome c Releasing Apoptosis Assay Kit
9. ApoSENSOR™ ATP Determination Kit & ADP/ATP Ratio Assay Kit
10. ApoGSH™ Glutathione Detection Kit
11. Nitric Oxide Fluorometric & Colorimetric Detection Kits
12. Annexin V-FITC, -Cy3, -Cy5, -EGFP, -PE, -Biotin Kits and Bulk Reagents
13. CaspSCREEN™ Flow Cytometric Caspase Assay Kit
14. CaspGLOW™ Active Caspase Staining Kits
15. CaspSELECT™ Caspase-3 & -7 Immunoassay Kits
16. Caspase-1,-2,-3,-5,-6,-8,-9,-10 Activity Assay Kits
17. Recombinant Human, Mouse, Rat Active Caspases & Procaspsases
18. Ready-to-use Caspase Inhibitors & Sets
19. Ready-to-use Caspase Substrates and Sets
20. Quick Apoptotic DNA Ladder Detection Kit
21. TUNEL-based DNA Fragmentation Assay Kits
22. Ready-to-use Apoptosis Inducers and Set
23. Quick Cell Proliferation Assay Kit
24. LDH-Cytotoxicity Assay Kit
25. Live-Dead Cell Staining Kit
26. Anti-Caspase-1,-2,-3,-4,-5,-6,-7,-8,-9,-10,-11,-12,-13,-14 & Active Caspases
27. Anti-Apaf-1, Bcl-2, Bcl-x, Bax, Bad, Bid, Bim/Bod, Mcl-1, Cytochrome c
28. Anti-DFF40, DFF45, FADD, p53, IRAK, Wee1, CD40, DcR2, cFLIP, CTGF
29. Many Other Apoptosis & Cell Signaling Related Antibodies
30. Anti-GFP, His, GST, and Other Tag Antibodies
31. Growth Factors, Cytokines & Chemokines (many)
32. Link-FAST™ 5 Minutes DNA Ligation Kit
32. Gel-FAST™ 20 Minutes Gel Staining/Destaining Kit
33. Protease Inhibitor Cocktail