

Recombinant Human Procaspase-9

CATALOG #: 1089P-5

AMOUNT: 5 µg

PURITY: ≥90% by SDS-PAGE

FORM: 2 µg/µl in 1 mM acidic acid

STORAGE:

Store at -70°C. Aliquot product into individual tubes and immediately refreeze at -70°C.
AVOID FREEZE/THAW CYCLES!

BACKGROUND DESCRIPTION:

Caspase-9 is a member of the caspase family of cysteine proteases. Caspase-9 exists in cells as an inactive proenzyme. Proteolytic processing of this inactive proenzyme generates the large and small subunits which, when assembled as a tetramer (a pair of heterodimers), form the active caspase. Cascades of caspase activation have been shown to be the key signal-transducing events in apoptosis. Caspase-9 is considered to be an upstream enzyme that initiates caspase cascade leading to apoptosis.

The recombinant Procaspase-9 (1089P-5) is produced by expression of a human cDNA in *E. coli* and further purified by multiple chromatographic steps.

APPLICATIONS:

Investigating mechanisms of the proteolytic processing and activation of procaspase-9.

NOTE: For research purpose only! Not to be used in human.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptotic Cell Isolation Kit

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay Kit
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Quality Antibodies for Apoptosis and Signal Transduction